

KRUUSE Livestock Equipment and Consumables

2' Edition





INDEX

Diagnostics	01 · 4
Gloves - Clothing	02 · 14
Injection	03 · 34
Identification	04 · 48
Wound Care - Bandages	05 · 50
Lubricants	06 · 54
Transfusion-Blood Collection	07 · 56
Hoof Treatment	08 · 58
Calf Rearing-Lactation	09 · 65
Magnets -Stomach Tubes & Trocars	10 · 68
Instruments-Desinfection	11 · 72
KRUUSE Obstetrics	12 · 84
Dehorning - Restraining	13 · 92
Feed Supplements	14 · 96
Cattle A.I. Misc.	15 · 98
Pig A.I.	16 · 104
Sheep Breeding Equipment	17 · 119



KRUUSE Digital Thermometer DIGI-Vet SC 1

- Temperature range: 32°C - 43.9°C
- Measuring accuracy: 0.1°C between 35.5°C to 42.0°
- Latest measuring result stored
- Turns off automatically after use
- Long temperature sensor for large animals

Cat. No	Description
291102	KRUUSE Digital Thermometer DIGI-Vet SC 12



KRUUSE DIGI-TEMP Express Thermometer

- Fast measuring: 10-15 sec.
- Temperature measurement between 32-42 °C
- Accuracy: +/- 0.10 °C
- Long temperature sensor for large animals

Cat. No	Description
291114	KRUUSE DIGI-TEMP Express Thermometer



KRUUSE Digital Thermometer, watertight

- Temperature range: 32°C - 43°C
- Latest measuring result stored
- Acoustic signal
- More than 200 hours of continuous operation on one battery 1.5 V

Cat. No	Description
291103	KRUUSE Digital Thermometer, watertight



KRUUSE Digital Thermometer w/flexible tip, waterproof

- Temperature range: 32°C - 43.9°C
- Measuring accuracy: 0.1°C between 35.5°C to 42.0°
- Latest measuring result stored
- Acoustic signal
- Turns off automatically after use

Cat. No	Description
291111	KRUUSE Digital Thermometer w/flexible tip, waterproof



Laboratory Thermometer

Temperature range -10/+110°C. Length 26 cm.

Cat. No	Description
291125	Laboratory thermometer



KRUUSE Thermometer Covers, 1000/pk

- For easy handling and protection of the thermometer
- Fits many standard digital thermometers

Cat. No	Description
299020	KRUUSE Thermometer Covers, 1000/pk



Littmann Stethoscope

A luxury model stethoscope. Known for its great quality and sound. Floating diaphragm enables listening to high and low frequency sounds by changing the pressure on the diaphragm.

Cat. No	Description
291204	Littmann Master Classic II for large animals. Blue tube 81 cm



Sprague Stethoscope

Dual head with separate double tubing. Supplied complete with 3 spare bells, 2 spare sets of ear pieces and 2 spare diaphragms. Specially designed for veterinary use, with 60 cm long tubing.

Cat. No	Description
291288	Stethoscope, Sprague black



Panascope Stethoscopes

Cat. No	Description
291260	Panascope economy model
291270	Panascope standard model



Thinklabs One

Features:

- The smallest, most powerful stethoscope in the world
- Fits in the palm of your hand
- Amplifies more than 100X
- Uses audio headphones

Description:

The icon of medicine has changed. Thinklabs One has a standard 3.5mm jack so you can connect high quality headphones instead of old tubing style binaurals. You get better sound, more comfort, a more compact stethoscope and freedom to use the headphones of your choice. The One includes studio quality In Ear Headphones, with excellent bass for heart sounds and crystal-clear treble for music enjoyment. Designed for high power output, the One has dedicated headphone amplifiers, one for each channel. Voltage doubling charge pumps provide extra headroom for low distortion reproduction of S1 and S2 heartbeat peaks. The headphone amplifiers have sufficient power to drive low-impedance audio headphones, producing 45.2mWatts. Competitive stethoscopes operate on lower voltage, use a single amplifier and high impedance speaker, resulting in clipped heartbeats at significantly lower power output. The One allow you to listen with multiple filter choices. You can also set the two filter choices you most often use, as "favorites" that are super quick to select. So you'll have both choice and speed. Ranging from very low frequency (low pitched sounds such as S3), to higher frequencies (lung sounds) and a few in between. The One uses a rechargeable Lithium Ion battery so there's no need to replace batteries. The USB-style charger cable is included, along with compact world-wide AC charger (100-240V). The battery capacity has been designed to require only one or perhaps two charge cycles per week. Thinklink provides a connection between One and your mobile or laptop devices to capture and send sounds via iPhone, iPad, Android devices and Mac/PC. Easily switch between listening to patients or playback recordings. Thinklink is included with One.

Includes:

- Thinklabs One Digital Stethoscope
- Thinklabs Studio Performance In-Ear headphones. (The One can be used with any high quality headphones.)
- Thinklink connection system for iPhone, iPad, Android, PC/Mac
- USB-type Charger 100-240V
- Pocket carrying case
- 2 Year Thinklabs Warranty

Cat. No	Description
291320	Thinklabs One



KRUUSE Urine Strips

Inexpensive and easy-to-read test. Coloured indicator squares on the strip can be compared with the colour chart on the label. Results appear in 60 seconds. Parameters: pH, nitrites, protein, glucose, ketones, bilirubin, urobilinogen, blood, leukocytes, specific gravity (For specific gravity we advise to use our refractometers).

Cat. No	Description
290798	KRUUSE Urine Strips 100/pk



Ketostix Ames

Urine stix, measures ketosis in urine. 50/pk
Measuring ketosis in the urine is more accurate than measuring ketosis in milk.

Cat. No	Description
290837	Ketostix ames 2880 50/pk



BHB milk

The test strip shows the concentration of BHB in milk.

BHB (Beta Hydroxy Butyrate) is a substance that in higher concentration is an indicator for ketosis.

Cat. No	Description
290839	BHB milk, 25/pk



KRUUSE BoDia Quick

Calf diarrhea poses a major health problem in herds, and is therefore a significant economic risk in cattle breeding. The causes of diarrhea are multiple. Infectious and non-infectious factors can trigger the onset of calf diarrhea. Non-infectious factors that may trigger diarrhea include, among others, hygienic conditions, the vitality of the calf's immune system and the physical condition of the calf.

Infectious factors include viruses, bacteria, single-cell parasites and other causative organisms. Some of the most common causative organisms are:

- Rotavirus
- Coronavirus
- Escherichia Coli K99
- Cryptosporidium parvum


These organisms, among others things, alter the structure of epithelial cells in the gastrointestinal tract, which interfere with normal digestion and absorption of nutrients and cause diarrhea. A calf with diarrhea excretes an increased volume of faeces with elevated content of water, causing disequilibrium in fluid balance. When a calf suffers from diarrhea, immediate treatment is required. It is imperative that the cause of the illness quickly and accurately is identified.

KRUUSE BoDia Quick

- Rapid diagnostics and screening on site
- Direct identification of 4 infectious pathogens
- Read the result after only 10 minutes
- User-friendly test-kit
- Storage at room temperature

Test procedure for liquid faeces

- 1) Pick up a faeces with the pipette
- 2) Open the purple lid of the sample test tube and remove the sampler
- 3) Add 3 drops of faeces into the test tube. If the faecal sample is very watery add up till 8 drops totally
- 4) Return the sampler to the tube and close it tightly. Shake the tube well and break off the pin in the opposite end of the test tube. Apply 3-4 drops of sample fluids into each of the sample wells by slightly squeezing the tube



Test procedure for grainy faeces or faeces with larger particles

- 1) Remove the purple lid of the test tube and use the sample stick, within in tube, for sampling
- 2) Stab the pin into three different places of the faeces
- 3) Return the sampler to the tube and close it tightly. Shake the tube well The clear liquid should change to a slightly brown/yellow colour
- 4) Break off the pin in the opposite end of the test tube. Apply 3-4 drops of sample fluids into each of the sample wells by slightly squeezing the tube



Cat. No	Description
296058	KRUUSE BoDia Quick 5/pk

KRUUSE PiDia Light

Rotavirus and Clostridium perfringens in pigs

Diarrhea can cause severe losses in livestock farming, especially during the first 3 weeks of life. Diarrhea generally leads to deterioration in individual animal health and can cause serious economic losses in herds.

Dominant pathogens vary between different herds. Diarrhea is usually caused by mixed infections. A thorough diagnosis is therefore of particular importance.

KRUUSE PiDia Light can reliably detect two of the most common pathogens causing diarrhea: Rotavirus and Clostridium perfringens. Problems in livestock management, such as husbandry, feeding and climate, can also cause or exacerbate diarrhea. In good livestock management, regular screening is recommended to check for possible increased concentrations of pathogens in order to identify and handle problems as fast and effectively as possible.

For example, if a high concentration of Clostridium perfringens (>10⁴ Clostridium perfringens/gram) is detected when screening, immediate preventive measures should be taken.


KRUUSE PiDia Light has been developed as an on-site diagnostic test for mixed infections causing diarrhea. Furthermore KRUUSE PiDia Light is also a screening instrument in herd management.

KRUUSE PiDia Light

- Rapid diagnostics and screening on site
- Read the result after 5 to 10 minutes
- Userfriendly test-kit
- Storage at room temperature

Test procedure for liquid faeces

- 1) Pick up a faeces with the pipette
- 2) Open the purple lid of the sample test tube and remove the sampler
- 3) Add 3 drops of faeces into the test tube. If the faecal sample is very watery add up till 8 drops totally
- 4) Return the sampler to the tube and close it tightly. Shake the tube well and break off the pin in the opposite end of the test tube. Apply 3-4 drops of sample fluids into each of the sample wells by slightly squeezing the tube



Test procedure for grainy faeces or faeces with larger particles

- 1) Remove the purple lid of the test tube and use the sample stick, within in tube, for sampling
- 2) Stab the pin into three different places of the faeces
- 3) Return the sampler to the tube and close it tightly. Shake the tube well The clear liquid should change to a slightly brown/yellow colour
- 4) Break off the pin in the opposite end of the test tube. Apply 3-4 drops of sample fluids into each of the sample wells by slightly squeezing the tube



Cat. No	Description
296059	KRUUSE PiDia Light, 5/pk



KRUUSE Müller Hinton Agar

Used for staphylococci, streptococci, micrococci and coliforms.
Often used for intestinal bacteria in large animal practice.
Can in many ways replace blood agar plates
Long shelf-life (versus blood agar)

Cat. No	Description
290614	Müller Hinton agar, 6/pk

KRUUSE Mannitol Salt Agar

For detecting staphylococci. Staphylococci are salt tolerant and can grow in up to 10% salt.

Contains red phenol as an indicator.
Pathogenic staphylococci: yellow zone around the colonies
Non-pathogenic staphylococci: small colonies without yellow zone



Cat. No	Description
290616	Mannitol salt agar, 6/pk

Glutavac Test

For the diagnosis of acute inflammation (increased fibrinogen concentration) and chronic inflammatory infections (increased immunoglobulin concentration) in large animal practice.
The GLUTAVAC test reagent is supplied in ready-to-use vacu tubes (10 ml) into which whole blood has to be drawn (like a standard blood sample).



Cat. No	Description
290785	GLUTAVAC test 10/pk



CMT Test Paddle

A four-well plastic paddle. Place milk from each quarter into each section. Add test liquid, agitate gently and evaluate. Marks indicate 2 and 3 ml.

Cat. No	Description
170365	CMT test paddle

Squeeze the handle and new test liquid comes out automatically.

Cat. No	Description
170390	CMT test paddle, X-spurt multitest, blue
170392	CMT test paddle, X-spurt multitest, white



CMT Test Liquid

Test liquid for the detection of mastitis. To be used with CMT test paddle. The test liquid indicates the cell count in milk.

Cat. No	Description
170366a	CMT liquid 1 l
170367	CMT liquid 5 l
170374	CMT liquid 10 l



KRUUSE Colostrum Densimeter

For quick evaluation of IgG content in bovine colostrum. A quick way of measuring the quality of the colostrum is by using the KRUUSE colostrum densimeter. The densimeter measures the specific gravity of the colostrum.

Cat. No	Description
290855	KRUUSE colostrum densimeter complete with measuring cylinder
290856	Spare densimeter
290857	Spare measuring cylinder, 250 ml



CMT Test Kit

An accurate and inexpensive method of determining the cell count in milk.

The kit contains:

- 500 ml dispenser bottle with dosing pump which gives exactly 3 ml test liquid
- CMT paddle
- Book with detailed instructions and pages for registration of test results

Test liquid is supplied separately.

Cat. No	Description
170362a	CMT test kit, English
170362	CMT test kit, Danish



290015

KRUUSE Refractometers

KRUUSE Colostrum Refractometer will very rapidly give a reliable colostrum measurement result, thus an indication if the colostrum quality is good enough to feed the newborn. Colostrum quality is measured in percentages with a cut-off value of 22%. Test the colostrum before feeding it to the newborn and use only good quality colostrum.

Cat. No	Description
290015	Colostrum refractometer, digital
290017	Colostrum refraktometer, Brix scale



Inoculating Sterile Needle

Disposable inoculating needle with eye, 10 x 10/pk.
Suitable for inoculation on agar plates etc.

Cat. No	Description
170370	Inoculating Sterile needle with eye 10 x 10/pk.



Thermocult incubator

Handy and inexpensive incubator for multi-purpose laboratory use - for example for warming up glassware and instruments. With thermometer. Temperature range 35-45 degrees Celcius, 230 V, 50 Hz.

Cat. No	Description
290370	Thermocult incubator



Vaginal Plexi Speculum

w/Xenonlight
(batteries not included)

Cat. No	Description
291000	30 x 250 mm w/Xenonlight
291001	40 x 350 mm w/Xenonlight
290002	50 x 420 mm w/Xenonlight



Polansky Speculum for Vaginal Examinations

The speculum has 3 blades. Has butterfly screws with 2 blades on top for wide exposure.

Nickel plated

Cat. No	Description
290990	Vaginal speculum Polansky 27 cm
290995	Vaginal speculum Polansky 35 cm



Metal Speculum for Vaginal Examinations

Nickel plated

Cat. No	Description
290985	Vaginal speculum for horses and cows, 27 cm metal duck bill type
290980	Vaginal speculum for horses and cows, 35 cm metal duck bill type
290979	Vaginal speculum for sheep 20 cm



Rondo Measuring Tape

For measuring slaughtering weight of pigs and cattle. Automatic recoil back into plastic case by pressing the button.
Length 2.5 m

Cat. No	Description
240589	Rondo combi weight measuring tape for pigs and cattle




Gloves

You have the choice between 3 different materials when you choose examination gloves: latex, vinyl or nitrile.

If you are allergic to latex we recommend nitrile gloves. If economy is your most important parameter, choose vinyl.

Nitrile gloves have become very popular as a milking glove in dairy farms due to its superior strength. Powder-free vinyl gloves are very popular at boar stations when collecting semen.



Measuring chart for Examination gloves

Instructions:

Size is indicated by the width of your hand just below your first knuckles. Place your right hand flat (palm down) on the guide with your fingers together and thumb spread. Your index finger should be flush against the dotted line, but not covering it.

**FOR MORE INFORMATION,
PLEASE SEE OUR SEPARATE GLOVE CATALOGUE.
ASK YOUR KRUUSE REPRESENTATIVE FOR YOUR FREE
COPY OR DOWNLOAD THE GLOVE CATALOGUE
AT WWW.KRUUSE.COM**



KRUTEX Latex Examination Gloves - Powder Free

High quality, latex examination gloves. Smooth, elastic and with high sensitivity. Non-sterile, ambidextrous. Powder free. Thickness at the fingers: 0.16 mm. Thickness at the palm: 0.15 mm. AQL: 1.5.

Features and Benefits:

- Easy dry hand donning and removal
- Powder free for extra allergen protection
- Palm textured for good wet and dry grip properties
- Extra thin for increased sensitivity and feel
- Free of Thiuams and Thiazoles
- Cuff strengthened by a rolled cuff
- Biodegradable environmentally friendly
- Hot water leached
- Single chlorination process
- Soft yet strong, well above minimum requirements
- Easy size recognition on box
- Colour coded for easy recognition

Cat. No	Description
260758	KRUTEX examination gloves, latex, PF, S, 100/pk
260759	KRUTEX examination gloves, latex, PF, M, 100/pk
260760	KRUTEX examination gloves, latex, PF, L, 100/pk
260769	KRUTEX examination gloves, latex, PF, XL, 100/pk



KRUTEX Latex Examination Gloves - Powdered

High quality, latex examination gloves. Smooth, elastic and with high sensitivity. Non-sterile, ambidextrous. Powder: Cornstarch. Thickness at the fingers: 0.16 mm. Thickness at the palm: 0.15 mm. AQL: 1.5.

Features and Benefits:

- Lightly powdered for easy donning and removal
- Smooth surface
- Extra thin for increased sensitivity and feel
- Free of Thiuams and Thiazoles
- Cuff strengthened by a rolled cuff
- Biodegradable environmentally friendly
- Hot water leached
- Easy size recognition on box
- Colour coded for easy recognition

Cat. No	Description
260761	KRUTEX examination gloves, latex, P, S 100/pk
260762	KRUTEX examination gloves, latex, P, M 100/pk
260763	KRUTEX examination gloves, latex, P, L 100/pk
260764	KRUTEX examination gloves, latex, P, XL 100/pk



KRUTEX Vinyl Examination Gloves - Powder Free

High quality vinyl examination gloves. Non-sterile. Non-spermicidal. For all-round examination purposes. Ambidextrous. Thickness at the fingers: 0.15 mm. Thickness at the palm: 0.12 mm. AQL: 1.5.

Features and Benefits:

- Easy donning and removal
- Cuff strengthened by a rolled cuff
- Extra thin for increased sensitivity and feel
- Super stretchy (fatigue free)
- Free of latex
- Free of protein
- Powder free for extra allergen protection
- Super soft providing superior user comfort
- Easy size recognition on box
- Colour coded for easy recognition

Cat. No	Description
260787	KRUTEX examination gloves, vinyl, PF, S, 100/pk
260788	KRUTEX examination gloves, vinyl, PF, M, 100/pk
260789	KRUTEX examination gloves, vinyl, PF, L, 100/pk
260790	KRUTEX examination gloves, vinyl, PF, XL, 100/pk



KRUTEX Vinyl Examination Gloves - Powdered

High quality vinyl examination gloves. Non-sterile. Non-spermicidal. For all-round examination purposes. Ambidextrous. Powder: Cornstarch. Thickness at the fingers: 0.15 mm. Thickness at the palm: 0.12 mm. AQL: 1.5.

Features and Benefits:

- Easy dry hand donning and removal
- Cuff strengthened by a rolled cuff
- Extra thin for increased sensitivity and feel
- Super stretchy (fatigue free)
- Free of latex
- Free of protein
- Super soft providing superior user comfort
- Easy size recognition on box
- Colour coded for easy recognition

Cat. No	Description
260791	KRUTEX examination gloves, vinyl, P, S, 100/pk
260792	KRUTEX examination gloves, vinyl, P, M, 100/pk
260793	KRUTEX examination gloves, vinyl, P, L, 100/pk
260794	KRUTEX examination gloves, vinyl, P, XL, 100/pk

KRUTEX Nitrile Clinic Gloves - Powder Free



A strong and durable glove. Non-sterile. Provides excellent protection against the most common drugs and chemicals used in animal production. Thickness at the fingers: 0.15 mm. Thickness at the palm: 0.13 mm. AQL: 1.5

Features and Benefits:

- Easy donning and removal
- Cuff strengthened by a rolled cuff
- Micro textured providing good wet and dry grip properties
- Extra thin for increased sensitivity and feel
- Good viral barrier properties
- Latex like stretch and comfort
- Softer than standard nitrile and vinyl
- Free of latex
- Free of powder
- Free of protein
- Easy size recognition on box
- Colour coded for easy recognition

Cat. No	Description
261131	KRUTEX nitrile examination gloves, powder free, XS, 100/pk
261132	KRUTEX nitrile examination gloves, powder free, S, 100/pk
261133	KRUTEX nitrile examination gloves, powder free, M, 100/pk
261134	KRUTEX nitrile examination gloves, powder free, L, 100/pk
261135	KRUTEX nitrile examination gloves, powder free, XL, 100/pk

KRUTEX Milking Gloves - Powder Free



Extra strong nitrile examination glove. Designed to be used as a milking glove. Can also be used as a clinic glove.

Features and Benefits:

- Easy gloving and removal
- Cuff strengthened by a rolled cuff
- Extra thin for increased sensitivity and feel
- Finger Textured Surface
- Free of Latex
- Free of Protein
- Free of Powder
- Good fit, extra smooth, silky feel
- Latex like stretch and comfort
- Good viral barrier properties
- Polymer coated powder-free
- Easy size recognition on box
- Colour coded for easy recognition

Cat. No	Description
261136	KRUTEX milking gloves, nitrile, PF, extra strong, S, 100/pk
261137	KRUTEX milking gloves, nitrile, PF, extra strong, M, 100/pk
261138	KRUTEX milking gloves, nitrile, PF, extra strong, L, 100/pk
261139	KRUTEX milking gloves, nitrile, PF, extra strong, XL, 100/pk
261141	KRUTEX milking gloves, nitrile, PF, extra strong, XXL, 100 /pk



KRUUSE Sleeve Cover

Light blue sleeve cover, 40x20 cm, one size. Light weight cover made of PE. Ideal in cases where liquids are used.

Cat. No	Description
260926	KRUUSE Sleeve Cover, 100/pk



Metal Stand for Disposable Gloves

Cat. No	Description
320232	Metal stand for gloves-holds 3 boxes
320234	Metal stand for gloves-holds 2 boxes






Surgical Gloves

Surgical gloves are in high demand in veterinary clinics and hospitals. Comfort, feel and protection all need to be excellent.

KRUUSE has for many years supplied the veterinary profession with KRUTEX surgical gloves, which in many countries are the vets' preferred choice. Together with leading glove experts we have upgraded the quality of KRUTEX surgical gloves.



Measuring Chart for Surgical and Examination Gloves

Instructions:

Size is indicated by the width of your hand just below your first knuckles. Place your right hand flat (palm down) on the guide with your fingers together and thumb spread. Your index finger should be flush against the dotted line, but not covering it.





KRUTEX Latex Surgical Gloves - Powder Free

The polymer coating ensures extra soft user comfort during long processes. The cuff is strengthened to provide excellent viral barrier properties while roll and tearing of the glove is prevented. The gloves are designed for single use, are powder-free and have been sterilised using gamma irradiation. The gloves are easy to put on and take off. Tested for chemical residue and irritation. Tested for cytotoxicity. Protein level: 50 ug/gram material. Biodegradable and eco-friendly.

Features and Benefits:

- Powder-free for extra allergic protection
- Good dry donning capabilities
- Micro-textured providing good grip for efficient instrument handling
- Hand specific with curved finger design for maximum dexterity
- Free of Thiurams and Thiazoles
- Low protein and undetectable chemical due to extensive leaching during manufacturing process
- Optimum thickness distribution for excellent tactile sensitivity
- Biodegradable environmentally friendly
- Excellent viral barrier properties
- Packed in premium PET/PE film pouch for added protection

Cat. No	Description
260958	KRUTEX surgical gloves, PF, sterile size 5.5 50 pairs
260961	KRUTEX surgical gloves, PF, sterile size 6.0, 50 pairs
260962	KRUTEX surgical gloves, PF, sterile size 6.5, 50 pairs
260963	KRUTEX surgical gloves, PF, sterile size 7.0, 50 pairs
260964	KRUTEX surgical gloves, PF, sterile size 7.5, 50 pairs
260965	KRUTEX surgical gloves, PF, sterile size 8.0, 50 pairs
260966	KRUTEX surgical gloves, PF, sterile size 8.5, 50 pairs
260967	KRUTEX surgical gloves, PF, sterile size 9.0, 50 pairs



KRUTEX Latex Surgical Gloves - Powdered



The polymer coating ensures extra soft user comfort during long processes. The cuff is strengthened to provide excellent viral barrier properties while roll and tearing of the glove is prevented. The gloves are designed for single use and have been sterilised using gamma irradiation.

- The gloves are easy to put on and take off
- Tested for chemical residue and irritation
- Tested for cytotoxicity
- Protein level: 50 ug/gram material
- Biodegradable and eco-friendly

Features and Benefits:

- Lightly powdered for excellent dry donning capabilities
- Micro-textured providing good grip for efficient instrument handling
- Hand specific with curved finger design for maximum dexterity
- Free of Thiurams and Thiazoles
- Low Protein and un-detectable chemical due to extensive leaching during manufacturing process
- Optimum thickness distribution for excellent tactile sensitivity
- Biodegradable environmentally friendly
- Excellent viral barrier properties
- Packed in premium PET/PE film pouch for added protection.

Cat. No	Description
260971	KRUTEX surgical gloves, P, sterile size 6.0, 50 pairs
260972	KRUTEX surgical gloves, P, sterile size 6.5, 50 pairs
260973	KRUTEX surgical gloves, P, sterile size 7.0, 50 pairs
260974	KRUTEX surgical gloves, P, sterile size 7.5, 50 pairs
260975	KRUTEX surgical gloves, P, sterile size 8.0, 50 pairs
260976	KRUTEX surgical gloves, P, sterile size 8.5, 50 pairs





VALUELINE Surgical Gloves - Powdered

Gloves from our economical line. Easy to put on and take off. Irritation tested.

Features and Benefits:

- Lightly powdered for excellent dry donning capabilities
- Hand specific with curved finger design for maximum dexterity
- Free of Thiurams and Thiazoles
- Low protein and undetectable chemical due to extensive leaching during manufacturing process
- Biodegradable environmentally friendly
- Microtextured providing good grip for efficient instrument handling

Cat. No	Description
260981	VALUELINE surgical gloves, P, sterile size 6.0, 50 pairs
260982	VALUELINE surgical gloves, P, sterile size 6.5, 50 pairs
260983	VALUELINE surgical gloves, P, sterile size 7.0, 50 pairs
260984	VALUELINE surgical gloves, P, sterile size 7.5, 50 pairs
260985	VALUELINE surgical gloves, P, sterile size 8.0, 50 pairs
260986	VALUELINE surgical gloves, P, sterile size 8.5, 50 pairs





KRUTEX Arm Length Gloves

KRUTEX Arm Length Gloves are used by veterinarians all over the world and are renowned for their extreme strength and very high sensibility. We have gloves available in 5 different qualities depending on the sensibility needed. The lower sensibility you need, the cheaper the glove.

In addition to the traditional glove design we also offer a neck strap version which has become increasingly popular. Not only will the neck strap hold the glove in place, it also provides added shoulder, armpit and chest protection for increased hygiene.

All products are manufactured according to ISO 527-3 (film tensile test).



KRUTEX Super Sensitive Examination Gloves

Our highest quality and sensibility.

Disposable arm length gloves, smooth, extremely strong and with ultra high tactile sensitivity. For rectal examinations, artificial insemination and obstetrics.

Cut diagonally for a better fit. Thickness 25 µm, 100% polyethylene.

Cat. No	Description
260693	KRUTEX super sensitive gloves, small, 90 cm, 100/pk
260700	KRUTEX super sensitive gloves, large, 90 cm, 100/pk
260687	KRUTEX super sensitive gloves, extra long 97 cm, 100/pk



KRUTEX Super Sensitive Examination Gloves with Shoulder Protector

Our highest quality and sensibility. Disposable arm length gloves with neck strap to provide shoulder and armpit coverage. Smooth, extremely strong and with ultra high tactile sensitivity. For rectal examinations, artificial insemination and obstetrics. Length 90 cm, thickness 25 µm, 100% polyethylene.

Cat. No	Description
260726	KRUTEX super sensitive gloves, with shoulder protector, L, 50/pk
260727	KRUTEX super sensitive gloves, with shoulder protector, S, 50/pk





KRUTEX Sonosleeve Examination Gloves

Specially developed 6 fingered glove with room for a linear, endorectal transducer in the 6th finger - all in one glove. This glove is ideal for rectal examination with ultrasound. A strong and smooth examination glove with high tactile sensitivity. Length 90 cm, thickness 25 mμ, 100% polyethylene.
Available with or without shoulder protector. Patented.

Cat. No	Description
260706	KRUTEX SonoSleeve glove, L, 100/pk
260707	KRUTEX SonoSleeve glove with shoulder protector, L, 50/pk



KRUTEX Super Touch Examination Gloves

High quality gloves with superior tactile sensitivity. Cut diagonally for a better fit.
For rectal examinations, artificial insemination and obstetrics.
Length 95 cm, thickness 25 mμ, 100% polyethylene.

Cat. No	Description
260711	KRUTEX super touch gloves, L, 100/pk



KRUTEX Super Touch Gloves with Shoulder Protector

High quality gloves with superior tactile sensitivity. Neck strap to keep in place to provide shoulder and armpit coverage for better hygiene.
For rectal examinations, artificial insemination and obstetrics.
Length 90 cm, thickness 25 mμ, 100% polyethylene.

Cat. No	Description
260728	KRUTEX super touch gloves with shoulder protector, S, 50/pk
260725	KRUTEX super touch gloves with shoulder protector, L, 50/pk



KRUTEX Soft Examination Gloves

Our most popular green disposable arm length glove. Soft and strong for multi purpose veterinary and AI use. Cut diagonally for a better fit. Length 95 cm, thickness 25 mμ, 100% polyethylene.

Cat. No	Description
260685	KRUTEX soft examination gloves, large, 100/pk



KRUTEX Soft Examination Gloves with Shoulder Protector

Our most popular green disposable arm length glove with neck strap to provide shoulder and armpit coverage. Soft and strong for multi purpose veterinary and AI use. Length 90 cm, thickness 25 mμ, 100% polyethylene.

Cat. No	Description
260724	KRUTEX soft examination gloves with shoulder protector, large, 50/pk



VALUELINE Examination Gloves

Our economy line glove. Long and durable examination glove. Generally used in cases where great sensibility is not required. Length 95 cm, thickness 25 mμ, 100% polyethylene.

Cat. No	Description
260690	VALUELINE examination glove, large, 100/pk





KRUTEX AGRO disposable examination gloves

Armlength examination gloves for rectal examinations, artificial insemination and obstetrics.
Thickness 20 mμ, Length 90 cm, 100% polyethylene.

Cat. No	Description
260718	KRUTEX AGRO disposable examination gloves, large, 100/pk



KRUTEX Obstetric Gown

Green, 69% PVC and 31% nylon polyamide.
Waterproof, with welded seams, lightweight quality gown, Velcro closing, rubber cuff ring in the sleeves. Comfortable to wear and individually adjustable. Size indicates length of the gown.

Cat. No	Description
260010	KRUTEX obstetric gown with Velcro, 120 cm
260020	KRUTEX obstetric gown with Velcro, 130 cm
260030	KRUTEX obstetric gown with Velcro, 140 cm
260040	KRUTEX obstetric gown with Velcro, 150 cm
260045	KRUTEX obstetric gown with Velcro, 160 cm



KRUTEX Obstetric Combi Suit

KRUTEX Obstetric Suit is shaped like a boiler suit. Green, 69% PVC and 31% nylon polyamide. Waterproof, with welded seams, lightweight quality, Velcro closing on the back and rubber cuff ring in the sleeves. Comfortable to wear.

Cat. No	Description
260050	KRUTEX obstetric combi suit, S, size 46
260051	KRUTEX obstetric combi suit, M, size 50
260052	KRUTEX obstetric combi suit, L, size 54
260053	KRUTEX obstetric combi suit, XL, size 58



KRUTEX Obstetric Jacket and Trousers

Green, 69% PVC and 31% nylon polyamide. Waterproof, with welded seams, lightweight quality. Velcro closing on the back and rubber cuff ring in the sleeves. Comfortable to wear.

Cat. No	Description
260055	KRUTEX obstetric jacket, XS
260056	KRUTEX obstetric jacket, S
260060	KRUTEX obstetric jacket, M
260070	KRUTEX obstetric jacket, L
260075	KRUTEX obstetric jacket, XL
260078	KRUTEX obstetric trousers, XS
260079	KRUTEX obstetric trousers, S
260080	KRUTEX obstetric trousers, M
260090	KRUTEX obstetric trousers, L
260095	KRUTEX obstetric trousers, XL



KRUTEX Obstetric Overalls

Green, 69% PVC and 31% nylon polyamide. Waterproof, with welded seams, lightweight quality. Comfortable to wear with adjustable braces.

Cat. No	Description
260096	KRUTEX obstetric overalls, S, size 46
260097	KRUTEX obstetric overalls, M, size 50
260098	KRUTEX obstetric overalls, L, size 54
260099	KRUTEX obstetric overalls, XL, size 58





KRUTEX Obstetric Apron

Green, 69% PVC and 31 % nylon polyamide.
Waterproof, with welded seams, lightweight quality,
tape closure in the neck.

Cat. No	Description
260450	KRUTEX obstetric apron 110 cm, green
260460	KRUTEX obstetric apron 120 cm, green
260470	KRUTEX obstetric apron 130 cm, green



KRUTEX Shoulder Protector

Green, 69% PVC and 31% nylon polyamide.
Waterproof, with welded seams, lightweight quality with rubber cuff ring.
Heavy duty and unisize.

Cat. No	Description
260560	KRUTEX shoulder protector, rubber sleeve plate, long version
260561	KRUTEX shoulder protector, rubber sleeve plate, short version



KRUTEX Obstetric Sleeve

Green, 69% PVC and 31% nylon polyamide.
Waterproof, with welded seams, lightweight quality with rubber cuff ring.
Heavy duty and unisize.

Cat. No	Description
260570	KRUTEX obstetric sleeve, short



KRUTEX Disposable Obstetric Gown

The Obstetric gown can be used for any obstetric procedure. A great protective gown.

5 selling points beavers

- Light green
- 0.03 mm thickness
- Waterproof
- Sleeves to be cut to desired length
- Available in 2 sizes

Cat. No	Description
260589	KRUTEX disposable obstetric gown, S, 25/pk
260609	KRUTEX disposable obstetric gown, L, 25/pk



KRUUSE Disposable Coveralls

Made of PP-spun bonded non-woven material with a hood. Ideal for farm visits.

Cat. No	Description
260608	KRUTEX disposable coveralls, green, S
260611	KRUTEX disposable coveralls, green, M
260612	KRUTEX disposable coveralls, green, L
260613	KRUTEX disposable coveralls, green, XL
260614	KRUTEX disposable coveralls, green, XXL
260607	KRUTEX disposable coveralls, green, XXXL



KRUTEX Disposable Polyethylene Gown, Green

The KRUTEX Disposable Gown is made of 100% polyethylene. Green, length 135 cm, with long sleeves.

Cat. No	Description
260615	KRUTEX disposable polyethylene gown, green, 20/pk

**KRUTEX Disposable Plastic Gown**

Green plastic gown with a tape closure in the back and at the wrists. Length 130 cm. This disposable coat offers superior protection and comfort.

Cat. No	Description
260610	KRUTEX disposable coat, green, non-sterile, 25/pk
260620	KRUTEX disposable coat, green, sterile, 25/pk

**KRUTEX Disposable Coat**

White plastic coat with press buttons in the front. Length 120 cm. This disposable coat offers superior protection and comfort.

Cat. No	Description
260605	KRUTEX disposable coat, 120 cm

**KRUTEX Disposable Apron**

White plastic apron made of 0.03 mm LD/LLD. Length 140 cm.

Cat. No	Description
260629	KRUTEX disposable apron, sterile, 25/pk



KRUTEX Disposable Apron

White plastic apron made of 0.03 mm LD/LLD.

Cat. No	Description
260627	KRUTEX disposable aprons with dispenser, 140 cm, 75/pk
260628	KRUTEX disposable aprons with dispenser, 160 cm, 75/pk
260626	Stainless steel rack for disposable aprons



260626



KRUTEX Disposable Apron, non-sterile

White plastic aprons made of 0.03 mm LD/LLD.
This disposable apron provides superior protection and comfort.

Cat. No	Description
260630	KRUTEX disposable aprons, 140 cm, 150/pk, non sterile



Ankle laces

Elasticated top

KRUTEX Disposable Boots

KRUTEX disposable boots are well-known due to high quality, strength and non-slip. No welding seams under the foot, comfortable fitting. 100% polyethylene.

Available in standard quality of 40 mμ (thickness) or in an extra heavy duty quality of 70 mμ (thickness).

Cat. No	Description
260637	Boot cover, long/thin with ankle laces, 25 pairs, L: 50 cm
260640	Boot cover, long/heavy with ankle laces, 10 pairs. L: 50 cm
260642	Boot cover, long/thick, elasticated top, 25 pairs, L: 53 cm



KRUTEX Non-slip Disposable Boot Cover

Boots with an elasticated top that stretches wide yet provides a snug fit to keep debris out. Heavy duty textured bottom that will reduce slipping. Boots are easy to put on. Length: 40 cm.

Cat. No	Description
260656	KRUTEX Non-slip disposable boot cover, XXL, 20/pk



KRUUSE Disposable Shoe Covers

Plastic shoe covers with elastic edge. Easy to put on. Breathable, soft, and lightweight. Full elastic band by ankles for a superior fit. Delivered in a box.

Cat. No	Description
260646	KRUUSE disposable shoe covers, blue, 100/pk




KRUTEX Practice Coats

Durable, for daily work in the clinic. Knee length, comfortable quality, 60% Dacron and 40% cotton, with snap fasteners. Non-iron. The size is equal to chest size.

Cat. No	Description
260115	KRUTEX practice coat, 100 cm, grey
260125	KRUTEX practice coat, 104 cm, grey
260135	KRUTEX practice coat, 108 cm, grey
260145	KRUTEX practice coat, 112 cm, grey
260155	KRUTEX practice coat, 116 cm, grey
260156	KRUTEX practice coat, 120 cm, grey





KRUUSE Disposable Needles

A range of sizes specially designed for the veterinary profession.
Extra sharp facet shaped needles with plastic hub, luer slip.
Paper-film peel pack.

Cat. No	Description	Colour
121288	Disp. Needle 0.4 x 20 mm, 27G x ¾", 100/pk	Grey
121289	Disp. Needle 0.45 x 25 mm, 26G x 1", 100/pk	Brown
121293	Disp. Needle, 0.5 x 10mm, 25Gx3/8", 100/pk	Orange
121270	Disp. Needle 0.5 x 16 mm, 25G x 5/8", 100/pk	Light Orange
121290	Disp. Needle 0.55 x 25 mm, 24G x 1", 100/pk	Purple
121271	Disp. Needle 0.6 x 16 mm, 23G x 5/8", 100/pk	Blue
121272	Disp. Needle 0.6 x 25 mm, 23G x 1", 100/pk	Light Blue
121286	Disp. Needle 0.6 x 30 mm, 23G x 1¼", 100/pk	Light Blue
121291	Disp. Needle 0.7 x 20 mm, 22G x ¾", 100/pk	Grey
121287	Disp. Needle 0.7 x 30 mm, 22G x 1¼", 100/pk	Grey
121275	Disp. Needle 0.8 x 10 mm, 21G x ¾", 100/pk	Light Green
121276	Disp. Needle 0.8 x 16 mm, 21G x 5/8", 100/pk	Light Green
121273	Disp. Needle 0.8 x 25 mm, 21G x 1", 100/pk	Light Green
121274	Disp. Needle 0.8 x 40 mm, 21G x 1½", 100/pk	Light Green
121292	Disp. Needle 0.9 x 13 mm,, 20G x ½", 100/pk	Yellow
121277	Disp. Needle 0.9 x 25 mm, 20G x 1", 100/pk	Yellow
121278	Disp. Needle 0.9 x 40 mm, 20G x 1½", 100/pk	Yellow
121281	Disp. Needle 1.1 x 16 mm, 19G x 5/8", 100/pk	White
121279	Disp. Needle 1.1 x 25 mm, 19G x 1", 100/pk	White
121280	Disp. Needle 1.1 x 40 mm, 19G x 1½", 100/pk	White
121282	Disp. Needle 1.2 x 40 mm, 18G x 1½", 100/pk	Pink
121283	Disp. Needle 1.25 x 25 mm, 18G x 1", 100/pk	Pink
121284	Disp. Needle 1.6 x 25 mm, 16G x 1", 100/pk	White
121285	Disp. Needle 1.6 x 40 mm, 16G x 1½", 100/pk	White



KRUUSE-Vet Needles

Reusable needles, made of stainless steel, sterilisable.
With extra strong luer lock hub for increased safety.



Cat. No	Description
120768	Vet Needles 0.7 x 4mm 22G x 1/8" luer lock, 12/pk
120607	Vet needles 0.8 x 10 mm 21G x 3/8", luer lock, 12/pk
120766	Vet Needles 0.9 x 6 mm 20G x 1/4" luer lock, 12/pk
120602	Vet needles 0.9 x 13 mm 20G x 1/2", luer lock, 12/pk
120827	Vet Needles 1.0 x 6 mm 19G x 1/4" luer lock, 12/pk
120603	Vet needles 1.0 x 10 mm 19G x 3/8", luer lock, 12/pk
120613	Vet needles 1.0 x 15 mm 19G x 5/8", luer lock, 12/pk
120614	Vet needles 1.0 x 20 mm 19G x 3/4", luer lock, 12/pk
120615	Vet needles 1.0 x 25 mm 19G x 1", luer lock, 12/pk
120612	Vet needles 1.0 x 30 mm 19G x 1 1/4", luer lock, 12/pk
120829	Vet Needles 1.2 x 6 mm 18G x 1/4" luer lock, 12 pk
120767	Vet Needles 1.2 x 10 mm 18G x 3/8" luer lock, 12/pk
120765	Vet Needles 1.2 x 13mm 18G x 1/2" luer lock, 12/pk
120722	Vet needles 1.2 x 15 mm 18G x 5/8", luer lock, 12/pk
120624	Vet needles 1.2 x 20 mm 18G x 3/4", luer lock, 12/pk
120723	Vet needles 1.2 x 25 mm 18G x 1", luer lock, 12/pk
120828	Vet Needles 1.2 x 30 mm 18G x 1 1/4" luer lock, 12/pk
120625	Vet needles 1.2 x 35 mm 18G x 1 3/8", luer lock, 12/pk
120622	Vet needles 1.2 x 40 mm 18G x 1 1/2", luer lock, 12/pk
120735	Vet needles 1.4 x 20 mm 17G x 1/2", luer lock, 12/pk
120728	Vet needles 1.4 x 25 mm 17G x 1", luer lock, 12/pk
120729	Vet needles 1.4 x 35 mm 17G x 1 3/8", luer lock, 12/pk
120732	Vet needles 1.4 x 40 mm 17G x 1 1/2", luer lock, 12/pk
120733	Vet needles 1.4 x 50 mm 17G x 2", luer lock, 12/pk
120743	Vet needles 1.6 x 13 mm 16G x 1/2", luer lock, 12/pk
120736	Vet needles 1.6 x 20 mm 16G x 3/4", luer lock, 12/pk
120744	Vet needles 1.6 x 25 mm 16G x 1", luer lock, 12/pk
120742	Vet needles 1.6 x 30 mm 16G x 1 1/4", luer lock, 12/pk
120745	Vet needles 1.6 x 40 mm 16G x 1 1/2", luer lock, 12/pk
120752	Vet needles 1.8 x 25 mm 15G x 1", luer lock, 12/pk
120753	Vet needles 1.8 x 40 mm 15G x 1 1/2", luer lock, 12/pk
120755	Vet needles 2.0 x 13 mm 14G x 1/2", luer lock, 12/pk
120759	Vet needles 2.0 x 20 mm 14G x 3/4", luer lock, 12/pk
120763	Vet needles 2.0 x 30 mm 14G x 1 1/4", luer lock, 12/pk
120824	Vet needles 2.0 x 40 mm 14G x 1 1/2", luer lock, 12/pk
120826	Vet needles 2.0 x 50 mm 14G x 2", luer lock, 12/pk




KRUUSE Disposable Needle with Aluminium Hub

Complete range of high quality injection needles. A range of sizes specially designed for the veterinary profession. Pyrogen-free and sterilized with ethylene oxide.

Cat. No	Description
121379	Disp. Needle with aluminium hub 1.27x25 mm, 18G x 1", 100/pk
121380	Disp. Needle with aluminium hub 1.6 x 25 mm, 16G x 1", 100/pk
121381	Disp. Needle with aluminium hub 1.6 x 38 mm, 16G x 1½", 100/pk
121384	Disp. Needle with aluminium hub, 2.0 x 25 mm, 14G x 1", 100/pk
121382	Disp. Needle with aluminium hub 2.0 x 38 mm, 14G x 1½", 100/pk
121383	Disp. Needle with aluminium hub 2.0 x 50 mm, 14G x 2", 100/pk

BOVIVET Disposable Needles

Disposable needles with extra strong nylon hub for heavy duty use.



Cat. No	Description
121412	Disp. Needle 0.8 x 80 mm, 21G x 3¼", luer lock 100/pk
121310	Disp. Needle 1.25 x 38 mm, 18G x 1½" luer lock 100/pk
121325	Disp. Needle 1.5 x 50 mm, 17G x 2" luer lock 100/pk
121330	Disp. Needle 1.6 x 25 mm, 16G x 1" luer lock 100/pk
121340	Disp. Needle 1.6 x 38 mm, 16G x 1½" luer lock 100/pk
121350	Disp. Needle 2.1 x 38 mm, 14G x 1½" luer lock 100/pk
121360	Disp. Needle 2.1 x 60 mm, 14G x 2½" luer lock 100/pk
121410	Disp. Needle 2.1 x 80 mm, 14G x 3¼" luer lock 100/pk

Safety Container for Sharp Objects

For used blades, needles and other sharp objects.



Cat. No	Description
122540	Container for used needles 0.5 l
122545	Container for used needles 2 l



KRUUSE Disposable Syringes


KRUUSE disposable syringes are available in 2 and 3 component versions and with centric or eccentric luer slip and luer lock attachment. We also offer a range of disposable syringes with silicone O-ring as well as disposable syringes with attached needles for easy use.

KRUUSE Disposable Syringe

Complete range of high quality disposable syringes, 2- and 3-component.

- The disposable syringes are available as centric and eccentric
- Sterilised by ethylene oxide
- Paper-film peel pack
- Extended graduation for added value
- Luer slip adaptor

Cat. No	Description
112380	Disp. Syringe centric 2-comp. 2-2,4 ml, 100/pk
112381	Disp. Syringe centric 2-comp. 5-5,8 ml, 100/pk
112382	Disp. Syringe eccentric 2-comp. 10-11 ml, 100/pk
112383	Disp. Syringe eccentric 2-comp. 20-23 ml, 50/pk
112372	Disp. Syringe centric 3-comp. 1 ml, 100/pk
112373	Disp. Syringe centric 3-comp. 2-3 ml, 100/pk
112374	Disp. Syringe centric 3-comp. 5-6 ml, 100/pk
112375	Disp. Syringe eccentric 3-comp. 10-12 ml, 100/pk
112376	Disp. Syringe eccentric 3-comp. 20-24 ml, 50/pk
112377	Disp. Syringe eccentric 3-comp. 30-34 ml, 50/pk
112378	Disp. Syringe eccentric 3-comp. 50-60 ml, 25/pk






VALUELINE Disposable Syringes

The range of VALUELINE syringes is an economic alternative to the other KRUUSE syringes.

- Sterilised by ethylene oxide
- Luer slip adaptor

Cat. No	Description
112442	Disp. Syringe centric 3-comp. 1 ml, 100/pk
112443	Disp. Syringe centric 3-comp. 2 ml, 100/pk
112444	Disp. Syringe centric 3-comp. 5 ml, 100/pk
112445	Disp. Syringe centric 3-comp. 10 ml, 100/pk
112446	Disp. Syringe eccentric 3-comp. 20 ml, 100/pk




KRUUSE Disposable Syringes Luer Lock

KRUUSE offers a complete range of high quality luer lock disposable syringes, 3-component.

- Sterilised by ethylene oxide
- Paper-film peel pack
- Extended graduation for added value
- Luer lock adaptor
- Centric tip

Cat. No	Description
112384	Disp. Syringe 3-comp. luer lock 2-3 ml, 100/pk
112385	Disp. Syringe 3-comp. luer lock 5-6 ml, 100/pk
112386	Disp. Syringe 3-comp. luer lock 10-12 ml, 100/pk
112387	Disp. Syringe 3-comp. luer lock 20-24 ml, 50/pk
112388	Disp. Syringe 3-comp. luer lock 30-34 ml, 50/pk
112389	Disp. Syringe 3-comp. luer lock 50-60 ml, 25/pk



EXTENDED GRADUATION

KRUUSE Disposable Syringes with Silicone O-ring

The KRUUSE Syringes with silicone O-ring are designed with a ring made of medical grade silicone rubber.


The cylinder piston is comprised and the thin silicone ring is mounted separately and is acting as a sliding seal.

The new design minimizes the contact area between the silicone and the solution which significantly reduces particulate contamination.

- Sterilised by ethylene oxide
- Paper-film peel pack
- Extended graduation for added value
- Luer slip adapter
- Centric tip

Cat. No	Description
112432	Disp. Syringe with silicone O-ring centric 3 comp. 1 ml 100/pk
112433	Disp. Syringe with silicone O-ring centric 3 comp. 2-3 ml 100/pk
112434	Disp. Syringe with silicone O-ring centric 3 comp. 5-6 ml 100/pk
112435	Disp. Syringe with silicone O-ring centric 3 comp. 10-12 ml 100/pk
112436	Disp. Syringe with silicone O-ring eccentric 3 comp. 20-24 ml 100/pk
112437	Disp. Syringe with silicone O-ring eccentric 3 comp. 50-60 ml 50/pk





KRUUSE Disposable Syringes with Fixed Needle

3-component syringes with fixed needle for easy use, with luer slip or luer lock adaptor.

- Extra sharp facet shaped needles
- With plastic hub, luer lock
- Paper-film peel pack
- Extended graduation for added value

Cat. No	Description	Colour
112416	Insulin syringe, 3-comp., 0.5 ml, 29G x 1/2", 100/pk	White
112417	Insulin syringe 3-comp., 1 ml, 29G x 1/2", 100/pk	White
112418	TBC syringe, 3-comp., 1 ml, 25G x 5/8", 100/pk	Orange
112419	Syringe 3-comp., 1 ml, 25G x 5/8", 100/pk	Orange
112420	Syringe 3-comp., 1 ml, 26G x 11/2", 100/pk	Brown
112438	Syringe 3-comp., 2-3 ml, 21G x 1/2", LL, 100/pk	Light Green
112430	Syringe 3-comp., 2-3 ml, 21G x 5/8", 100/pk	Light Green
112423	Syringe 3-comp., 2-3 ml, 21G x 11/2", 100/pk	Light Green
112431	Syringe 3-comp., 2-3 ml, 22G x 3/4", 100/pk	Dark Grey
112422	Syringe 3-comp., 2-3 ml, 22G x 11/4", 100/pk	Dark Grey
112421	Syringe 3-comp., 2-3 ml, 23G x 1", 100/pk	Blue
112447	Syringe 3-comp., 5-6 ml 18G x 11/2" 100/pk	Pink
112425	Syringe 3-comp., 5-6 ml, 21G x 11/2", 100/pk	Light Green
112424	Syringe 3-comp., 5-6 ml, 22G x 11/4", 100/pk	Dark Grey
112427	Syringe 3-comp., 10-12 ml, 20G x 11/2", 100/pk	Yellow
112426	Syringe 3-comp., 10-12 ml, 21G x 11/2", 100/pk	Light Green
112428	Syringe 3-comp., 10-12 ml, 22G x 11/2", 100/pk	Dark Grey
112429	Syringe 3-comp., 20-24 ml, 18G x 11/2", 100/pk	Pink
112439	Syringe luer lock, 3-comp., 1 ml, 25G x 5/8", 100/pk	Orange
112441	Syringe luer lock, 3-comp., 2-3 ml, 21G x 11/2", 100/pk	Light Green
112448	Disp. syringe with needle, 3 component, 2/3 ml luer slip, 21G x 5/8", 100/pk	Light Green
112430	Syringe luer lock, 3-comp., 2-3 ml, 22G x 3/4", 100/pk	Dark Grey
112440	Syringe luer lock, 3-comp., 2-3 ml, 22G x 11/4", 100/pk	Dark Grey

Closing Caps for Syringes

Cat. No	Description
112395	Closing caps for syringes 100/pk, non-sterile
112396	Closing caps for syringes 100/pk, sterile







KRUUSE Reusable Syringes

Our range of high quality reusable syringes. Careful cleaning and maintenance will ensure a long life span.



KRUUSE Plexi Syringe

Unbreakable, high quality syringe with TPX plexiglass material. Sterilisable up to 130°C. Exchangeable metal luer lock nozzle.

Cat. No	Description
112318	Plexi syringe, luer lock, 5 ml
112280	Plexi syringe, luer lock, 10 ml
112285	Plexi syringe, luer lock, 20 ml
112286	Plexi syringe, luer lock, 30 ml
112290	Plexi syringe, luer lock, 50 ml
112281	Plexi syringe, luer lock, 10 ml, with graduated piston rod
112287	Plexi syringe, luer lock, 20 ml, with graduated piston rod
112295	Plexi syringe, luer lock, 30 ml, with graduated piston rod
112289	Plexi syringe, luer lock, 50 ml, with graduated piston rod
112321	Washer for Plexi syringe 5 ml
112291	Washer for Plexi syringe 10 ml
112292	Washer for Plexi syringe 20 ml
112298	Washer for Plexi syringe 30 ml
112293	Washer for Plexi syringe 50 ml
112326	Silicone oil 18 ml for plexi syringe
112325	KRUUSE Oral Cannula f/ Plexi Syringe



BOVIVET PLUS Plexi Revolver Syringe

Semi-automatic revolver syringe, 50 ml, with unbreakable plexi-glass barrel and luer lock

- Made of polyamide and fiberglass
- Weighs only 280 g, ergonomic design
- Supplied in a strong plastic case
- Dosing: 1-2-3-4-5 ml

Cat. No	Description
112650	BOVIVET PLUS plexi revolver syringe luer lock 50 ml



HSW ECO-MATIC® Automatic Multi-dose Syringe

Durable high strength plastic syringe. Perfect ergonomically designed. Accurate dosage. Can be disinfected in boiling water without loss of performance or accuracy. Each Syring is supported with sparepart kit

Cat. No	Description
113121	HSW ECO-MATIC® automatic syringe with bottle attachment, 0.3 ml, luer lock
113138	Eco-Matic automatic syringe 1 ml, LL, w/bottle attachment
113122	HSW ECO-MATIC® automatic syringe with bottle attachment 2 ml, luer lock
113125	HSW ECO-MATIC® automatic syringe with bottle attachment 5 ml, luer lock
113139	Eco-Matic automatic syringe 1 ml, LL, w/tube attachment
113112	ECO-MATIC automatic syringe 2ml, LL incl. tubing
113114	O-ring for HSW ECO-MATIC® 2 ml
113115	ECO-MATIC automatic syringe 5ml, LL incl. tubing
113120	O-ring for HSW ECO-MATIC® 5 ml



HSW ECO-MATIC® with Bottle Mount Universal Collars

Color coded bottle collars for easy reference and identification of the right adapters for the different sized bottles Easy use, durable design, safe and secured bottle fitting For use with all medication bottles of different sizes and shapes (10, 50, 100, 200, 250 or 500 ml and larger) depending on the collar used: white 20 mm, blue 30 mm, green 33 mm

Cat. No	Description
113181	HSW ECO-MATIC® 1 ml Luer Lock with bottle attachment for thread bottles
113182	HSW ECO-MATIC® 2 ml Luer Lock with bottle attachment for thread bottles
113183	HSW ECO-MATIC® 5 ml Luer Lock with bottle attachment for thread bottles
113184	Set of Draw-off Caps Ø 20, 30, 33 mm white, blue, green



HSW ECO-MATIC® TWIN – Automatic self-filling syringe with two separate barrels

Premium-quality plastic syringe with durable plastic barrels. Ergonomic, leight weight design for ease of use. Simultaneous vaccination of two products independently. Bottle mount and/or tubing connection. Adjustable barrel orientation for comfort. Easy adjustable barrel position. Various fixed dose barrel combinations available

Cat. No	Description
113104	HSW ECO-MATIC TWIN® 0,5ml/0,5ml
113106	HSW ECO-MATIC TWIN® 1ml/1ml
113107	HSW ECO-MATIC TWIN® 2ml/2ml
113108	HSW ECO-MATIC TWIN® 1ml/2ml
113109	HSW ECO-MATIC TWIN® 0,5ml/2ml



HSW ECO-MATIC® Draw-off Caps

Cap for use in draw-off with the Eco-Matic automatic syringe. Attractive design with a protective cover on the tip and a connector with two diameters for connection of several tube types. Can also be used for most syringes with tubing attachment.

Cat. No	Description
113116	Draw-off cap, HSW ECO-MATIC®, white, 20 mm for 50ml and 100ml bottles 50/pk
113117	Draw-off Cap, HSW ECO-MATIC®, blue, 30 mm for 250ml and 500ml bottles /50pk
113118	Draw-off Cap, HSW ECO-MATIC®, green, 33 mm for 500ml bottles or bigger 50/pk



HSW MULTI-MATIC® Syringes

Unbreakable multi-dose semi-automatic syringe. Accurate and easy dose setting. Low weight and ergonomically designed. Autoclavable at 134°C/273°F.

Cat. No	Description
112840	HSW MULTI-MATIC® revolver syringe, luer lock, 25 ml
112841	HSW MULTI-MATIC® revolver syringe, luer lock, 50 ml
112850	Cylinder 25 ml LL for HSW MULTI-MATIC® Syringe
112851	Cylinder 50 ml LL for HSW MULTI-MATIC® Syringe
112852	O-ring for 25 ml HSW MULTI-MATIC® Syringe
112853	O-ring for 50 ml HSW MULTI-MATIC® Syringe

HSW ROUX-REVOLVER® Syringes

Semi-automatic multidose syringe. Heavy duty with metal body, rubber plunger, glass barrel, luer lock.



Cat. No	Description
112615	HSW ROUX-REVOLVER® syringe luer lock 10 ml
112635	HSW ROUX-REVOLVER® syringe luer lock 30 ml
112655	HSW ROUX-REVOLVER® syringe luer lock 50 ml
112661	HSW ROUX-REVOLVER® syringe 50 ml with extension piece (35 cm)
112730	Spare barrel for 10 ml HSW ROUX-REVOLVER® syringe
112745	Spare barrel for 30 ml HSW ROUX-REVOLVER® syringe
112740	Spare barrel for 50 ml HSW ROUX-REVOLVER® syringe
112768	Washers and piston for 10 ml HSW ROUX-REVOLVER® syringe
112715	Washers and piston for 30/50 ml HSW ROUX-REVOLVER® syringe



HSW Henke-Ject® TBC Syringe

Luer lock 1.8 ml, 0.1 ml graduation. Especially developed for the requirements of tuberculosis examination. Exact and accurate dosing – 0.1ml per lever trigger.

Cat. No	Description
113605	HSW Henke-Ject® TBC syringe luer lock 1.8 ml, 0.1 ml graduation



Socorex Automatic Self-filling Syringes

The Socorex® self-refilling syringes with spring-loaded plunger and three-way valve system are designed for serial precision injections. Simple conception, perfect balance and ergonomic design provide ease of use, operator comfort and effortless fieldwork.

Cat. No	Description
110674	Socorex aut. syringe luer lock 0.5 ml with bottle attachment
110675	Socorex aut. syringe luer lock 1 ml with bottle attachment
110679	Socorex aut. syringe luer lock 2 ml with bottle attachment
110680	Socorex aut. syringe luer lock 5 ml with bottle attachment
110682	Socorex aut. syringe luer lock 10 ml with bottle attachment
110730	Socorex aut. syringe luer lock 1 ml with tubing attachment
110735	Socorex aut. syringe luer lock 2 ml with tubing attachment
110740	Socorex aut. syringe luer lock 5 ml with tubing attachment
110745	Socorex aut. syringe luer lock 10ml with tubing attachment
110751	Socorex Twin aut. syringe luer lock 0.5 ml with tubing attachment
110752	Socorex Twin aut. syringe luer lock 1 ml with tubing attachment
110753	Socorex Twin aut. syringe luer lock 2 ml with tubing attachment
110754	Socorex Twin aut. syringe luer lock 5 ml with tubing attachment





Socorex Accessories

Cat. No	Description
110672	Barrel and plunger for Socorex syringe 0.5 ml
110677	Barrel for Socorex syringe 1 ml
110678	Barrel for Socorex syringe 2 ml
110681	Barrel for Socorex syringe 5 ml
110683	Barrel for Socorex syringe 10 ml
110671	Set of washers for Socorex automatic syringe 0.5 ml
110676	Set of washers for Socorex syringe 1 ml
110684	Set of washers for Socorex aut. syringe 2 ml
110685	Set of washers for Socorex aut. syringe 5 ml
110708	Set of washers for Socorex aut. syringe 10 ml
110688	Piston rod for Socorex aut. syringe 0.5 ml and 1 ml
110689	Piston rod for Socorex aut. syringe 2 ml
110690	Piston rod for Socorex aut. syringe 5 ml
110687	Spring for Socorex aut. syringe 0.5 and 1 ml
110692	Spring for Socorex aut. syringe 2 ml
110693	Spring for Socorex aut. syringe 5 ml
110698	Piston for Socorex aut. syringe 1 ml
110699	Piston for Socorex aut. syringe 2 ml
110700	Piston for Socorex aut. syringe 5 ml
110747	Socorex feeding needle 19.5 cm
110670	Tubing with suction needle and sinker for Socorex automatic syringe
110712	Vial holder for Socorex automatic dosing syringe
110669	Vial holder complete with feeding needle and tube connector



KRUUSE Drenching Guns

For oral administration of liquids in ruminants.
Adjustable dosage - easy cleaning and maintenance
- accessories included.

Cat. No	Description
113880	KRUUSE drenching gun 30 ml non-automatic
113881	KRUUSE drenching gun 60 ml non-automatic
113884	KRUUSE drenching gun 200 ml non-automatic
113882	KRUUSE drenching gun 30 ml automatic
113883	KRUUSE drenching gun 60 ml automatic



BOVIVET Automatic Drenching Gun

Cat. No	Description
	Heavy duty drencher with metal body and plexi glass barrel. The dosing can be set to any volume. The drencher is delivered complete with a strong 5 l plastic container with shoulder straps and tubing.
113302	BOVIVET drenching gun automatic 30 ml including 5 l. water container
113303	Graduated glass barrel 30 ml
113304	Rear washer for drenching gun 30+50 ml
w112305	Silicone for BOVIVET plexi revolver syringe
112326	Silicone oil 18 ml for plexi revolver syringe
113308	HSW Drench-Matic automatic syringe 30 m
113309	Plastic container for HSW Drench-Matic automatic syringe, 2.5 l



HSW Drench-Matic Drenching Gun

Automatic drenching gun for administration of solutions and suspensions. For oral administration and pour-on. Light weight ergonomic aluminum handle and high strength plastic barrel. Easy handling and cleaning.

Cat. No	Description
112308	HSW Drench-Matic drenching gun without plastic bag pack
112309	HSW plastic bag pack, with funnel, carrying strap, closing cap and draw-off cap





Marking

Correct and clear identification and marking af livestock is important. KRUUSE offers several options for this purpose: Marking sticks, marking sprays, tatooing equipments.



BOVIVET Superior Marking Stick

Superior marking sticks made of special waxes, paraffin oil and high quality dyes. Temperature resistant and highly colour intensive. With screw bottom.

Cat. No	Description
240300	BOVIVET superior marking stick, blue, 10/pk
240301	BOVIVET superior marking stick, red, 10/pk
240302	BOVIVET superior marking stick, green, 10/pk
240303	BOVIVET superior marking stick, black, 10/pk
240304	BOVIVET superior marking stick, purple, 10/pk



KRUUSE Marking Stick

An economic alternative to the well-known BOVIVET Superior marking sticks. KRUUSE marking sticks contain approx. 20% more wax than average sticks on the market. The KRUUSE marking sticks do not contain animal fat. With screw bottom.

Cat. No	Description
240313	KRUUSE marking stick green 10/pk
240314	KRUUSE marking stick red 10/pk
240315	KRUUSE marking stick blue 10/pk
240316	KRUUSE marking stick black 10/pk



KRUUSE FDA Marking Spray

Quick and convenient marking of animals in all weather conditions. Dries instantly and has long durability. Spray at a distance of 30-40 cm by pressing the nozzle lightly. Approved according to Danish Food Safety Regulations.

Cat. No	Description
240393	KRUUSE FDA marking spray, red, 400 ml, 12/pk
240394	KRUUSE FDA marking spray, green, 400 ml, 12/pk
240395	KRUUSE FDA marking spray, blue, 400 ml, 12/pk

KRUUSE Marking Spray

Quick and convenient marking of animals in all weather conditions.
Dries instantly and has long durability.
Spray at a distance of 30-40 cm by pressing the nozzle lightly.



Cat. No	Description
240416	KRUUSE marking spray, red, 200 ml, 12/pk
240425	KRUUSE marking spray, green, 200 ml, 12/pk
240411	KRUUSE marking spray, blue, 200 ml, 12/pk
240415	KRUUSE marking spray, red, 400 ml, 12/pk
240420	KRUUSE marking spray, green, 400 ml, 12/pk
240410	KRUUSE marking spray, blue, 400 ml, 12/pk
240419	KRUUSE marking spray, red, 500 ml, 12/pk
240424	KRUUSE marking spray, green, 500 ml, 12/pk
240418	KRUUSE marking spray, blue, 500 ml, 12/pk



Tattoo Forceps, Single Row for Pigs

With double action forceps tattooing can be performed more accurately while using less force. Double action, 1 row of figures 1 x 5.

Cat. No	Description
240471	Tattoo forceps for one row of five 10 mm figures



Tattoo Figures

For 10 mm tattoo forceps cat No. 240471.

Cat. No	Description
240474	Alphabet A-Z, 10 mm
240475	Numbers 0-9, 10 mm



Tattoo Hammer

For easy tattooing of pigs

Cat. No	Description
240541	Tattoo hammer without figures
240542	Set of figures for tattoo hammer (240541)



KRUUSE Vet-Flex

High quality flexible cohesive bandage. Strong, comfortable and user-friendly. Easy to tear. The Vet-Flex bandage has minimal shrinkage. Each roll is wrapped individually. Dispenser box of 4 or 10 rolls.

Cat. No	Description
160728	KRUUSE Vet-Flex 2.5 cm x 4.5 m, red, 10/pk
160729	KRUUSE Vet-Flex 5 cm x 4.5 m, white, 10/pk
160730	KRUUSE Vet-Flex 5 cm x 4.5 m, red, 10/pk
160744	KRUUSE Vet-Flex 5 cm x 4.5 m, green, 10/pk
160745	KRUUSE Vet-Flex 5 cm x 4.5 m, blue, 10/pk
160746	KRUUSE Vet-Flex 5 cm x 4.5 m, black, 10/pk
160731	KRUUSE Vet-Flex 7.5 cm x 4.5 m, green, 10/pk
160747	KRUUSE Vet-Flex 7.5 cm x 4.5 m, red, 10/pk
160748	KRUUSE Vet-Flex 7.5 cm x 4.5 m, blue, 10/pk
160749	KRUUSE Vet-Flex 7.5 cm x 4.5 m, black, 10/pk
160732	KRUUSE Vet-Flex 7.5 cm x 4.5 m, white, 10/pk
160733	KRUUSE Vet-Flex 10 cm x 4.5 m, white, 10/pk
160734	KRUUSE Vet-Flex 10 cm x 4.5 m, red, 10/pk
160735	KRUUSE Vet-Flex 10 cm x 4.5 m, blue, 10/pk
160736	KRUUSE Vet-Flex 10 cm x 4.5 m, green, 10/pk
160737	KRUUSE Vet-Flex 10 cm x 4.5 m, yellow, 10/pk
160738	KRUUSE Vet-Flex 10 cm x 4.5 m, black, 10/pk
160739	KRUUSE Vet-Flex 10 cm x 4.5 m, ass. of neon colours, 10/pk
160742	KRUUSE Vet-Flex 15 cm x 15 m, red, 4/pk
160743	KRUUSE Vet-Flex 15 cm x 15 m, blue, 4/pk
160726	KRUUSE Vet-Flex 15 cm x 15 m, black, 4/pk

VALUELINE Cohesive Flex Bandage

Economy range of flexible, cohesive and comfortable bandages. Easy to tear, no scissors needed. Suitable for all species. An alternative to our KRUUSE Vet-Flex bandages.



Cat. No	Description
164010	VALUELINE cohesive flex bandage, 5 cm x 4.5 m, red, 10/pk
164000	VALUELINE cohesive flex bandage, 5 cm x 4.5m, white, 10/pk
164100	VALUELINE cohesive flex bandage 7.5 cm x 4.5 m, white, 10/pk
164120	VALUELINE cohesive flex bandage 7.5 cm x 4.5 m, green, 10/pk
164200	VALUELINE cohesive flex bandage 10 cm x 4.5 m, white, 10/pk
164210	VALUELINE cohesive flex bandage 10 cm x 4.5 m, red, 10/pk
164220	VALUELINE cohesive flex bandage 10 cm x 4.5 m, green, 10/pk
164230	VALUELINE cohesive flex bandage 10 cm x 4.5 m, blue, 10/pk
164240	VALUELINE cohesive flex bandage 10 cm x 4.5 m, black, 10/pk



KRUUSE Vet-Plast Elastic Adhesive Bandage

Adhesive bandage of white elastic cotton with a yellow application guide line in the center. Good ventilation. Reverse winding - adhesive outer side to facilitate application.

Can be used to fix a dressing in place and as support for sprains, joints and fractures.

Cat. No	Description
161821	KRUUSE Vet-Plast 5 cm x 4.5 m, 10/pk
161831	KRUUSE Vet-Plast 7.5 cm x 4.5 m, 10/pk
161836	KRUUSE Vet-Plast 10 cm x 4.5 m, 10/pk



KRUUSE Wound Plast Spray

Wound plaster spray. Transparent, water-resistant protective film without tar. Environmentally friendly. No CFC's.

Cat. No	Description
161010	KRUUSE Wound Plast spray 200 ml



KRUUSE Wound Plast Spray with Tar

Sterile wound plaster spray. Transparent, water-resistant protective film with tar. Fly-repelling. Environmentally friendly. No CFC's.

Cat. No	Description
161020	KRUUSE Wound Plast spray with Tar, 200 ml



KRUUSE Aluminium Wound Spray

An astringent and antibacterial spray, which promotes wound healing. Acts as a wound plaster, which allows the wound to breathe while keeping the wound clean.

For external use on animals only. Contains micronized aluminium. Apply to affected area 1-2 times daily.

Cat. No	Description
161035	KRUUSE Aluminium Wound spray 200 ml



KRUUSE Blue-Spray

A disinfectant spray for bite wounds, cuts and other wounds. Can be used for horses, dogs, cats, pigs and cattle.

Cat. No	Description
161060	KRUUSE Blue-Spray 250 ml with atomizer, 6/pk
161070	KRUUSE Blue-Spray refill 250 ml, 6/pk
161080	KRUUSE Blue-Spray, 200 ml aerosol



BOVIVET Ointment Spray

Zinc oxide ointment spray. Protects the skin against sweat, urine, water and other harmful substances. Environmentally friendly. No CFC's.

Cat. No	Description
161100	BOVIVET Ointment Spray, 200 ml, skin protection ointment



KRUSAN Zinc Ointment 200 ml

Water-repellent unscented protective ointment. Immediately after application, KRUSAN zinc ointment forms an oil film, which protects the skin against water, sweat, urine, faeces and other secretions.

KRUSAN zinc ointment is made as a thin zinc paste which is easy to distribute and is easily absorbed by the skin. Although the ointment is hydrophobic, it lets the skin breathe. KRUSAN zinc ointment has a low pH (4.5) to maintain the skin's natural acidity which is crucial to the skin's defence mechanism against ammonia formed by urine and faeces.

The active substance zinc oxide has a drying effect and relieves skin redness and irritation. This ointment also contains allantoin and chamomile to calm the skin.

Cat. No	Description
325076	KRUSAN zinc ointment, 20%, 200 ml



Tensoplast

Elastic adhesive bandage. Widely used for claw and hoof treatment.

Cat. No	Description
160820	Tensoplast 5 cm x 4.5 m
160830	Tensoplast 7.5 cm x 4.5 m
160835	Tensoplast 10 cm x 4.5 m



Omniplast 2.5 cm x 5 m

A hypoallergenic tape, adhering to the skin without risking allergic reactions.

Easy to tear off, and can be used to fix medical instruments and dressings on necrotic tissue.

Cat. No	Description
169006	Omniplast 2.5 cm x 5 m, 12/pk
169007	Omniplast 5 cm x 5 m, 6/pk



BOVIVET Teat bandage

Unique polyester bandage. Protects wound against dirt, easy to mold and apply, cohesive. Can be applied to wet teats. Allergy tested.

Cat. No	Description
160954	BOVIVET Teat bandage, 6 cm x 4.5 m



BOVIVET Gel

The reference brand in the industry, with protective film for lubrication of instruments, arms and hands for births, vaginal and rectal examinations. Ready for use.

Cat. No	Description
180520	BOVIVET gel 500 ml
180525	BOVIVET gel 1000 ml
180528	BOVIVET gel 5000 ml
180529	BOVIVET gel 10 l



VALUELINE Examination Gel

Economic alternative to the famous BOVIVET gel. Good quality lubricating gel for vaginal and rectal examinations. High viscosity. The tube facilitates the introduction of the gel into the vagina and the rectum. After use the bottle is closed by bending the tube onto a small knob on the screw cap.

Cat. No	Description
180545	VALUELINE examination gel 1000 ml
180548	VALUELINE examination gel 5000 ml



BOVIVET Gel Granulate

Mix with water to make a lubricating gel. It can be used for rectal/vaginal examination and insemination. 12 g of gel granulate is mixed with 1 l water.

Cat. No	Description
180511	BOVIVET gel granulate bottle, 350 g
180512	BOVIVET gel granulate refill, 350 g
340625	Steril gel





121559



121582



121634



121556



121568



121647



121574



121575

BD Vacutainers and Needles

Cat. No	Description
121559	BD Vacutainer 2 ml plain, 100/pk
121582	BD Vacutainer 4 ml plain, 100/pk
121634	BD Vacutainer serum clot activator 10 ml (367896), 100/pk
121556	BD Vacutainer 2 ml EDTA, 100/pk
121568	BD Vacutainer 6 ml EDTA, 100/pk
121647	BD Vacutainer 10 ml EDTA, 100/pk
121600	BD Vacutainer Heparin Hemogard w/sodium 5 ml, 100/pk
121587	BD vacutainer, natrium heparin, 10 ml, 100/pk
121574	BD Vacutainer needle 18Gx1 S/S, 100/pk
121575	BD Vacutainer needle 20Gx1 M/S, 100/pk
121525	BD Vacutainer needle holder, 25/pk



230604

BOVIVET Infusion

Disposable needles with nylon hub, luer lock

Cat. No	Description
121410	BOVIVET disposable needles 2.1 x 80 mm, luer lock 100/pk
230604	BOVIVET infusion set with vented spike for bovine use



121410



Infusor Combi

Handy and convenient disposable combi-set that allows blood transfusion or fluid therapy on large animals over a long period without constant supervision. Sterile packed.

Please note: Add 125 ml CPDA-1 to 875 ml blood to get 1000 ml for blood transfusion. This anticoagulant will give whole blood a shelf life of 35 days.

Cat. No	Description
230510	Infusor Combi set 5 l. 5/pk
230520	Spare Infusor plastic bottle 5 l, sterile, 8/pk
230483	CPDA-1 anti-coagulant 500 ml



Blood transfusion Set

Convenient disposable combi-set, allows blood transfusion or fluid therapy on large animals over a longer period without constant supervision. Sterile packed.



Cat. No	Description
230515	Infusor and irrigating set with needle 2.0 x 160 mm, 25/pk



KRUUSE Hoof Glue

Fast setting two-component polyurethane adhesive for cattle

KRUUSE Hoof Glue has been designed as a treatment for bovine hoof pathologies. It is used for bonding a wooden orthopedic block under the healthy hoof. Thus, the injured hoof is raised clear of the ground. That allows an immediate recovery of the mobility, a fast healing and the resumption of the normal milk yield in a short term.

Cat. No	Description
220600	KRUUSE Hoof Glue, 200 ml

Accessories for KRUUSE Hoof Glue

- It is highly recommended to use protective gloves ex. KRUTEX Nitrile Gloves

Cat. No	Description
220601	KRUUSE Hoof Glue mixing tips, 10/pk
220602	KRUUSE Hoof Glue application gun for 200 ml
220490	Demotec wooden blocks, 110 mm
220491	Demotec wooden blocks XL, 130 mm



Demotec 90

Convenient and inexpensive kit for hoof treatment on bovines. Simple to use. Fast reduction of lameness, no reduction of milk production or weight of animal.

Complete set for 2 treatments including: wooden blocks, powder, liquid, measuring cups and spatulas.

Cat. No	Description
220450	Demotec kit 2 sets
220455	Demotec kit 12 sets
220460	Demotec 90 liquid 500 ml
220470	Demotec powder, 1000

Demotec 95

Same excellent features as for DEMOTEC 90. But DEMOTEC 95 offers in addition a superior moulding consistency that does not stick to your hands and can be used immediately after being mixed. Supplied as a complete set for 2 treatments. Easy to use. Complete sets for 2 or 14 treatments including: wooden blocks, powder, liquid, measuring cups and spatulas.



Cat. No	Description
220492	Demotec 95 kit 2 sets
220493	Demotec 95 kit 14 sets
220494	Demotec 95 liquid 500 ml
220495	Demotec 95 powder, 1000 g

Demotec Easy Bloc

A new system for which the design is based on the anatomy of the claw. Minimizes lameness. Fast, safe and easy to use. Demotec Easy Bloc hoof shoes 4 or 12 sets (2 or 6 right and 2 or 6 left). Not available in UK.

Clean the healthy hoof removing dirt and loose horn. Pour the liquid into the hoof shoe and sprinkle powder onto the liquid. Mix the two components and then apply the hoof shoe. Pressure can be applied to the hoof in less than 10 minutes. The process can be accelerated by using a heat gun.



Cat. No	Description
220530	Demotec Easy Bloc hoof shoes 4 sets (2 right, 2 left)
220531	Demotec Easy Bloc hoof shoes 12 sets (6 right, 6 left)
220535	Demotec Easy Block XXL 12 units, 6 right, 6 left
220532	Demotec Easy Bloc powder 500 g
220533	Demotec Easy Bloc liquid 250 ml

Demotec Futura Pad

FuturaPad is unique. One surface of the wooden block employed is covered with a flexible reactive Pad, which overcomes any unevenness of the sole and maximises adhesion over the full surface area. The reactive Pad is shaped correctly according to anatomy and the function of the claw and is affixed to the wooden block in the optimal position. FuturaPad is reactive. Pad is activated by a special glue.



Cat. No	Description
220550	Demotec Futura Pad hoof shoes, 2/pk
220551	Demotec Futura Pad hoof shoes, 10/pk
220552	Demotec Futura Pad hoof shoes, left, 10/pk
220553	Demotec Futura Pad hoof shoes, right, 10/pk

Demotec Hoof Check

Used as control in connection with hoof trimming. Used to control the five following parameters: Controlling the angle between the sole and the anterior hoof wall, measuring the length of the anterior hoof wall, comparing the height of the inside and outside hooves, ensuring that optimum contact between the sole and the ground can be achieved and controlling the slope of the sole.



Cat. No	Description
220565	Demotec hoof check

Demotec Wooden Blocks

To be fastened with SureBond. To be glued on the healthy part of the claw for immobilization of the sick part of the claw.



Cat. No	Description
220490	Demotec wooden blocks, 110 mm, 150/pk
220491	Demotec wooden blocks XL, 130 mm, 125/pk



BOVIVET Hoof Bandage

BOVIVET hoof bandage is a complete bandage ready to apply. It is made from a water-repellent material that keeps liquids out. Its shape and elastic materials are designed to keep the stocking in place. A treatment pad placed inside the stocking makes dosing easy, such as salicylic acid, once the stocking is put in place. Applying additional doses is possible and the effect of the treatment can be checked during the process. The stocking remains in place, and you can be sure there will be no complications. The stocking is easy to remove and is available in 2 sizes: XL- fits most types of cows · L - fits heifers and small cow breeds.

Cat. No	Description
220411	BOVIVET Hoof Bandage, L, 10/pk
220412	BOVIVET Hoof Bandage XL, 10/pk



KRUUSE Hoof-Gel

Contains: Salicylic acid 38% and Aloe Vera
Prepare the hoof by washing it thoroughly with a soapy solution or use hydrogen peroxide and dry it. Apply the hoof gel to the hoof and apply a BOVIVET hoof bandage which will withstand water and moisture. The bandage must be removed after a maximum of 48 hours. In case of Equine Thrush the procedure can be repeated as required.

Cat. No	Description
220402	KRUUSE Hoof-Gel with Salicylic acid 38%, 500 ml



BOVIVET Hoof-Dressing

Dressing for the treatment of digital dermatitis and other hoof ailments. Used ideally in conjunction with e.g. Kruuse Hoof Gel with salicylic acid. Fixate the dressing on the hoof to help the preparation to work in a clean environment.


Cat. No	Description
220405	BOVIVET protective hoof-dressing, 50/pk



Hoof RaspsCat. No

For treatment of approximately 1000 cows. Reduces dust and heat generation.

Cat. No	Description
220064	Hoof rasp disc tungsten 115 mm, coarse
220065	Hoof rasp disc tungsten 115 mm, fine



Demotec Hoof Discs

This hoof milling cutter plate has been developed in cooperation with hoof specialists.

Cat. No	Description
220067	Demotec DL hoof trimming disc
220068	Knife for Demotec DL hoof trimming disc
220069	Mounting unit for Demotec hoof trimming disc



Hoof Knives

Broad-bladed stainless steel, length of blades 80 mm. Standard quality.

Cat. No	Description
220052	Hoof knife standard right
220053	Hoof knife standard left
220054	Hoof knife standard double edge



KRUUSE Hoof Knives

Stainless, very sharp, length of blades 80 mm. With hard wood handle.

Cat. No	Description
220035	KRUUSE hoof knife narrow bladed, right
220042	KRUUSE hoof knife broad bladed right
220043	KRUUSE hoof knife broad bladed left
220044	KRUUSE hoof knife broad bladed, double
220045	KRUUSE hoof knife, narrow bladed, left



Aesculap Hoof Knives

High quality for professional users.

Cat. No	Description
220007	Aesculap hoof knife right
220008	Aesculap hoof knife left



Swiss hoof knife

With wooden handle.

Cat. No	Description
220026	Swiss hoof knife



Sharpening Tools

For hoof knives and scrapers.

Cat. No	Description
220070	Sharpening stone standard
220075	Sharpening stone Hauptner
220077	Sharpening steel for hoof knife, pocket model



Hoof Trimmer

Double action, for easy use.

Cat. No	Description
220211	Hoof trimmer short, standard 41 cm
220216	Spare blades with screws (for 220211), 2/pk



Hoof Cutter

High quality for professional users. Double action.

Cat. No	Description
220230	Hoof cutter, 60 cm



Hoof Rasp Cat. No

Straight, 35 cm, extra deep and sharp rasp blow.

Cat. No	Description
220130	Hoof rasp straight 350 mm



Hoof Hammer

Cat. No	Description
220160	Hoof hammer with plastic head



Hoof Tar

Used in connection with hoof treatment. Makes the bandage water-repellant. In metal bucket with handle.

Cat. No	Description
220430	Hoof tar 2.5 L



Tar Bandage

5 cm x 25 m, waterproof, double adhesive, for hoof treatment.

Cat. No	Description
160990	Tar bandage 4.5 cm x 25 m



BOVIVET Nylon Hoof Bandage

Bandage for inter-digital cleft with Velcro fastening.

Cat. No	Description
220058	BOVIVET hoof bandage for inter-digital cleft



Vet-Vat Cotton

For use in claw treatment and other purposes.

Cat. No	Description
161110	Vet-Vat cotton wool 100 g





Drenchers, Feeding Tubes, Trocars, Magnets and Infusion


KRUUSE offers veterinarians in large animal practice well known and thoroughly tested products such as the BOVIVET Calf Drencher, BOVIVET Magnets, different types of trocars and several different feeding tubes as well as other products to help in the daily work with livestock.



230548 / 230554



230545 / 230553



230547 / 230601



230549 / 230603



230558 / 230559

BOVIVET Calf Drencher

The BOVIVET Calf Drencher is available with a rigid and a flexible probe and the container holds 2 or 4 litres.

The Calf Drencher is a quick and safe method to administer liquid therapy of calves and for electrolytic calf treatment and it can be operated by one person.

The large container is primarily used when initially administering 4 l colostrum to the calf.

Two-position valve for adjusting the flow of raw milk:

- If the valve is in position 1, the liquid runs slowly.
- If the valve is in position 2, the liquid runs faster.

The 2 litre container will empty on average in 1 minute and the 4 litre container will empty on average in 2 minutes.

Cat. No	Description
230545	BOVIVET calf drencher, 2 liter with rigid probe, 10/pk
230548	BOVIVET calf drencher, 2 litre, with flexible probe, 10/pk
230553	BOVIVET calf drencher, 4 litre, with rigid tube, 5/pk
230554	BOVIVET calf drencher, 4 litre, with flexible tube, 5/pk
230547	BOVIVET flexible tube for calf drencher, 2 liter
230601	BOVIVET flexible tube for calf drencher, 4 litre
230549	BOVIVET rigid tube for calf drencher, 2 litre
230602	BOVIVET rigid tube for calf drencher, 4 litre
230552	Bottle cap for calf drencher, 2 litre
230603	Bottle cap for calf drencher, 4 litre
230558	BOVIVET teat for calf drencher bottle, 2 litre
230559	BOVIVET teat for calf drencher bottle, 4 litre



Calf Fluid Feeder

With 3 l bag and flexible probe. The bag is made of non-toxic, heavy grade plastic. Tear-shaped ball is designed so that the oesophageal probe cannot be inserted into the trachea. Smooth plastic tube prevents needless irritation. The reusable feeder has convenient loop for hanging.

Cat. No	Description
230550	Calf Fluid Feeder with 3 litres bag and flexible probe



170823



170822



170826



170825



170835



170821

Plastic Buckets With Feeding Teats

Cat. No	Description
170823	Plastic bucket with 5 feeding teats
170825	Plastic bucket with 1 feeding teat
170822	Teat for plastic bucket (170823)
170821	Plastic valve for 170823 and 170825
170826	Extra feeding teat for 170825 and 170823, calves
170835	Bottle teat medium



KRUUSE Colostrum Densimeter

For quick evaluation of IgG content in bovine colostrum. A quick way of measuring the quality of the colostrum is by using the KRUUSE colostrum densimeter. The densimeter measures the specific gravity of the colostrum.

Cat. No	Description
290855	KRUUSE colostrum densimeter complete with measuring cylinder



290015



290017

KRUUSE Refractometers

KRUUSE Colostrum Refractometer will very rapidly give a reliable colostrum measurement result, thus an indication if the colostrum quality is good enough to feed the newborn. Colostrum quality is measured in percentages with a cut-off value of 22%. Test the colostrum before feeding it to the newborn and use only good quality colostrum.

Cat. No	Description
290015	Colostrum refractometer, digital
290017	Colostrum refraktometer, Brix scale





BOVIVET Ruminal Magnet

Our well-known BOVIVET ruminal magnets are recommended by veterinarians worldwide due to their excellent function. The magnets have the correct size, shape and weight to offer maximum protection against and treatment of traumatic indigestion in cattle. The correct length and deep grooves in the plastic cage provide extra protection from sharp, uneven pieces.

Cat. No	Description
142030	BOVIVET ruminal magnet, 10/pk
142035	BOVIMINI ruminal magnet, 10/pk



BOVIVET Ruminal Magnet, Extra Powerful

Treatment and prevention of traumatic reticulitis. Tests have showed significant reduction in occurrence of traumatic reticulitis in cows when ruminal magnets are used. The BOVIVET Mag II has a higher pulling resistance than a conventional magnet.

Cat. No	Description
142040	BOVIVET Mag II, extra powerful ruminal magnet 10/pk





KRUUSE Ceramic Magnets

Presented in a strong Styrofoam box with each magnet in a separate compartment for maximum protection.
Magnets are easy to pick up as they do not stick.

Cat. No	Description
142037	KRUUSE Ceramic magnets 60 x 20 x 15 mm, 12/pk



BOVIVET Magnet Introducer

For use with BOVIVET magnets for cows.

Cat. No	Description
142045	BOVIVET magnet introducer metal
142050	BOVIVET magnet introducer plastic, diameter 3 cm
142070	BOVIMINI magnet introducer metal



Metal Detector

The AF 103 is developed for security checks, but it is also eminently suitable for many other purposes, such as finding wall ties, pipes and cables in walls and other structures, scanning luggage, and especially for finding swallowed objects.
One control operation incorporating on/off and tuning 6" (15 cm) waterproof loop
Audio output with built-in speaker
Lightweight and strong mechanical construction
Pinpointing ability
Uses one PP3 9 volt battery
Detects all metals


Cat. No	Description
141753	Metal detector AF 103 incl. battery 9V



Thygesen Stomach Tube for Cows

Thygesen stomach tube, stainless steel. Screw-on-ends available, small and large cup, beet cutting loop.

Cat. No	Description
190010	Stomach tube Thygesen, 18 x 1800 mm
190050	Beet cutting loop for Thygesen tube
190030	Cup 34 mm for Thygesen tube
190020	Cup 31 mm for Thygesen tube



Thygesen Stomach Tube for Calves

Thygesen stomach tube, stainless steel. Screw-on-ends available, cup and beet cutting loop.

Cat. No	Description
190150	Stomach tube for calves, 12 x 1200 mm
190160	Cup for calf tube
190170	Beet cutting loop for calf tube



BOVIVET Stomach Tube

Complete with perforated plastic tubing, heavy duty nylon protecting tube, large and small cup, tubing for pump.

Cat. No	Description
190110	BOVIVET stomach tube, 15 x 4000 mm



BOVIVET Stomach Probang System Mark II

The tube is made from polyurethane. Polyurethane is stronger than PVC and the flexibility does not change by variation of temperatures. Polyurethane also has a very smooth surface and tolerates most of the common disinfectants.

Cat. No	Description
190180	BOVIVET stomach probang system Mark II, excl. pump
190181	Spare tube for BOVIVET probang system Mark II
190182	Large cup for BOVIVET rumen system MK2, 40 mm
190183	Small cup for BOVIVET rumen system MK2, 32 mm
190184	Protecting mouth piece for BOVIVET Rumen System MK2
190185	Quick coupler for BOVIVET rumen system MK2



Special Stomach Pump

Stomach pump for 190180. The pump is double-acting and nickel-plated. Capacity of approx. 20 l per minute. (700 ml per pump beat). The tube is connected directly to the pump by means of snap-action coupling. Total length 59 cm. Length of pump cylinder 40 cm and inner diameter of cylinder 4.77 cm.


Cat. No	Description
190190	Special Stomach Pump for BOVIVET Stomach Probang Mark II
190191	BOVIVET O-ring for stomach pump (190190)

190400



KRUUSE Trocar

Cat. No	Description
190430	KRUUSE Trocar 110 mm, plastic stylet
190400	KRUUSE Buff screw trocar 135 mm, metal stylet



Trocars

For fast relief of gas accumulated in the rumen.

Cat. No	Description
190365	Trocar 4 x 100 mm, stainless steel for sheep and calves with one cannula
190310	Trocar 7 x 120 mm, stainless, with one cannula and wooden handle
190395	Trocar 10 x 150 mm, stainless, with one cannula and wooden handle
190320	Spare cannula for 190310
190396	Spare cannula for 190395
190370	Trocar 3 x 135 mm stainless steel
190368	Trocar 4 x 120 mm with metal handle and cannula
190369	Spare cannula for 190368
190325	Trocar 7 x 120 mm, with metal handle and cannula
190326	Spare cannula for 190325



KRUUSE Sterile Scalpel Blades

Manufactured from the finest surgical and razor blade grade steel. A high percentage of carbon gives particularly good grinding and honing properties, facilitating the manufacture of high quality edges.

The higher carbon content also enables each blade to be heat treated to achieve the perfect balance between hardness and strength. The material also has good corrosion resistance properties and conforms to ISO 7153 and BS 5194 standards. KRUUSE scalpel blades fit all common brands of scalpel handles that conform to the International Standard ISO.

Cat. No	Description
141114	A KRUUSE sterile scalpel blade No 10, stainless steel, 100/pk
141115	B KRUUSE sterile scalpel blade No 11, stainless steel, 100/pk
141116	C KRUUSE sterile scalpel blade No 12, stainless steel, 100/pk
141117	D KRUUSE sterile scalpel blade No 15, stainless steel, 100/pk
141118	E KRUUSE sterile scalpel blade No 20, stainless steel, 100/pk
141119	F KRUUSE sterile scalpel blade No 21, stainless steel, 100/pk
141132	G KRUUSE sterile scalpel blade No 22, stainless steel, 100/pk
141141	H KRUUSE sterile scalpel blade No 23, stainless steel, 100/pk
141182	I KRUUSE sterile scalpel blade No 24, stainless steel, 100/pk

KRUUSE Disposable Scalpels

The scalpel handle is made from ABS plastic providing greater rigidity and body to the scalpel therefore ensuring a superior finish. All materials are latex free.

The scalpel blades are manufactured from the finest surgical and razor blade grade steel. A high percentage of carbon gives particularly good grinding and honing properties, facilitating the manufacture of high quality edges. The higher carbon content also enables each blade to be heat treated to achieve the perfect balance between hardness and strength. The material also has good corrosion resistance properties and conforms to ISO 7153 and BS 5194 standards.



Cat. No	Description
140963	KRUUSE disposable scalpel no 22, 10/pk
140964	KRUUSE disposable scalpel no 23, 10/pk
140965	KRUUSE disposable scalpel no 24, 10/pk



KRUUSE Scalpel Handles

Made of stainless steel.

Cat. No	Description
141010	KRUUSE Handle No 3
141110	KRUUSE Handle No 4



KRUUSE Carbon Steel Blades

Surgical scalpel blades made from carbon steel.

Carbon steel is somewhat harder than stainless steel. This may produce a very sharp scalpel blade. However, the edge is not as durable as in stainless steel.

Cat. No	Description
141100 A	KRUUSE carbon steel blades no 10, sterile, 100/pk
141101 B	KRUUSE carbon steel blades no 15, sterile, 100/pk
141102 C	KRUUSE carbon steel blades no 22, sterile, 100/pk
141103 D	KRUUSE carbon steel blades no 23, sterile, 100/pk
141104 E	KRUUSE carbon steel blades no 24, sterile, 100/pk

Caesarian Knife



Cat. No	Description
141170	Special knife for Caesarian, 50/pk

Autopsy Knives



Cat. No	Description
140844	Autopsy knife 20 cm pointed, butcher's model
140845	Autopsy knife 21.5 cm curved, skinning knife
140846	Autopsy knife 15 cm curved, skinning Knife
140847	Autopsy knife 15 cm pointed, boning Knife



Bleeding Knife

Cat. No	Description
140848	Knife 21 cm for bleeding large animals , 6/pk
140850	Sharpening steel 25 cm

KRUUSE Post Mortem Kit

- Consists of:
- 140840, Stainless chopper
 - 140845, Autopsy knife curved, 21.5 cm, stainless steel, with plastic handle
 - 140846, Skinning knife curved, 15 cm, stainless steel, with plastic handle
 - 140847, Boning knife, pointed, 15 cm, stainless steel, with plastic handle
 - 140850, Sharpening steel, 25 cm
 - 140855, Bone saw, stainless stell
 - 140856, Spine wrench, stainless steel
 - 140857, Brain knife, stainless steel
 - 140858, Cartilage knife
 - 140859, Bone chisel, stainless steel
 - 140860, Postmortem hammer
 - 131205, Bone scissors, stainless steel, 23 cm
 - 131135, Bowel scissors, 21 cm
 - 130430, Nursing scissors, 14 cm
 - 130455, Straight scissors B/P, 16 cm
 - 140240, Rochester pean, straight, 16 cm
 - 240588, Tape measure 2 m
 - 250170, Metal box with tray
 - 141164, Scalpel No 6
 - 141165, Scalpel No 9
 - 141166, Scalpel No 18

Cat. No	Description
140830	KRUUSE post mortem kit





140855



140650



140840



140856



140853

Post Mortem Tools

Cat. No	Description
140855	Bone saw, stainless steel, blade 14", approx. 35.5 cm
140650	KRUUSE wound retractor Weitlaner 14 cm
140840	Meat chopper
140856	Spine wrench
140857	Brain knife standard quality
140859	Bone chisel standard quality
140860	Post mortem hammer
140853	Blade for bone saw 14" model, 35.5 cm approx.

140857



140859



140860



Surgical Forceps

Stainless steel.


Cat. No	Description
140710	Surgical forceps 13 cm standard quality
140716	Surgical forceps 15 cm standard quality
140718	Surgical forceps 18 cm standard quality



Surgical Forceps

Stainless steel.


Cat. No	Description
140720	Anatomic forceps 13 cm standard quality
140725	Anatomic forceps curved 13 cm standard quality
140717	Anatomic forceps 15 cm standard quality
140721	Anatomic forceps 18 cm straight, standard quality



Straight Scissors

Stainless steel.

Cat. No	Description
130445	Scissors straight 14 cm pointed/blunt, standard quality
130447	Scissors straight 14 cm blunt/blunt, standard quality
130449	Scissors straight 14 cm pointed/pointed, standard quality
130455	Scissors straight 16 cm pointed/blunt, standard quality
130456	Scissors straight 16 cm blunt/blunt, standard quality
130457	Scissors straight 16 cm pointed/pointed, standard quality
130629	Scissors straight 18 cm blunt/pointed, standard quality
130630	Scissors straight 18 cm blunt/blunt standard quality
130631	Scissors straight 18 cm pointed/pointed, standard quality



Curved Scissors

Stainless steel.

Cat. No	Description
130020	Cooper scissors 14 cm, curved pointed/blunt, standard quality
130120	Cooper scissors 14 cm, curved blunt/blunt, standard quality
130215	Cooper scissors 14 cm, curved pointed/pointed, standard quality
130216	Cooper scissors 16 cm, curved pointed/blunt, standard quality



Suture Clips

Used for skin sutures instead of suture material.

Cat. No	Description
151905	Suture clips 7.5 mm, 100/pk
151909	Suture clips 11 mm, 100/pk
151910	Suture clips 12 mm, 100/pk
151920	Suture clips 14 mm, 100/pk
151930	Suture clips 16 mm, 100/pk
151940	Suture clips 18 mm, 100/pk
151950	Suture clips 20 mm, 100/pk
151960	Suture clips 22 mm, 100/pk
152020	Suture clips 29 mm, 100/pk



Applicator for Suture Clips

Clips forceps for skin agraffes. Aluminium, for 29 mm agraffes. MICHEL suture clips applicator. Also for removing suture clips, stainless steel.

Cat. No	Description
152010	Applicator for suture clips 16 cm
151820	Michel suture clips applicator, 12.5 cm



Disposable Staplers

Cat. No	Description
152127	Signet 35W stapler disposable
152129	Royal 35W stapler disposable
152141	Henry Schein disposable skin stapler, 35 regular (3.9 x 59 mm), 5/pk
152142	Henry Schein disposable skin stapler, 35 wide (4.9 x 7,2 mm), 5/pk



Skin Staple Removal Forceps

Cat. No	Description
152126	Skin staple removal forceps
152143	Henry Schein disposable staple remover



KRUUSE LDA Toggle Sutures

To be used for non-surgical fixation of the left displaced abomasum.

Cat. No	Description
152175	KRUUSE LDA toggle sutures, stainless, 2/pk





Grymer Sterner Suture Set

To be used for non-surgical fixation of the left displaced abomasum.

Cat. No	Description
152180	Suture set Grymer Sterner 2/pk
152185	Trocar with cannula for inserting Grymer Sterner toggle suture
152190	Pressure relieving plates for Grymer Sterner toggle suture, 6/pk



Castrators for Bloodless Castration

With cord stop and double action handles.

Cat. No	Description
200151	Bloodless castrator 23 cm standard quality, Burdizzo type
200131	Bloodless castrator 30 cm standard quality, Burdizzo type
200136	Bloodless castrator 40 cm standard quality, Burdizzo type
200141	Bloodless castrator 48 cm standard quality, Burdizzo type



Schick Injector Razor

Schick Injector razor, only for large animals.
Used for pre-operation preparation.

Cat. No	Description
142133	Razor for large animals
142134	Razor blades, 10/pk



BLITZ Slaughter Pistol

BLITZ slaughter pistol. For fast and safe destruction of animals.
Cartridges in various strengths (see colour indications)

Cat. No	Description
240210	BLITZ slaughter pistol




BLITZ Cartridges

Cartridges for BLITZ slaughter pistol, 50 pcs. Red, for very large bulls. Green, for pigs and calves. Yellow, for horses and cows. Blue, for bulls.

Cat. No	Description
240215	Cartridges for Blitz, red 50/pk
240220	Cartridges for Blitz, green 50/pk
240230	Cartridges for Blitz, yellow 50/pk
240240	Cartridges for Blitz, blue 50/pk



240243



240250



240255

BLITZ Spare Parts

Cat. No	Description
240243	Bolt for Blitz pistol 205
240248	Cleaning kit for Blitz
240250	Blitz set for breech, spring No 215, rubber socket No 216
240254	Head and release button for Blitz slaughter pistol
240255	Blitz complete striker set



240248



KRUSAN Liquid Hand Soap

KRUSAN liquid hand soap is a mild and gentle soap to be used when frequent hand wash is required. Special attention to nails and roots of the nails are required in hand wash.

Cat. No	Description
325010	KRUSAN liquid hand soap, 1000 ml for wall dispenser
325040	KRUSAN liquid hand soap, 300 ml bottle including pump
325043	KRUSAN hand soap, 5 l
325063	KRUSAN shampoo for normal hair, 5 l



KRUSAN Hand Disinfecting Gel

KRUSAN surgical hand disinfection gel with 85% alcohol is effective against microorganisms. The gel formulation ensures the optimum duration of action in the hard to reach areas round nails and roots of the nails. Glycerine added as skin care.

Cat. No	Description
325020	KRUSAN hand disinfecting gel, 85%, 1000 ml for wall dispenser



KRUSAN Wall Dispenser

Robust dispenser in stainless steel

Cat. No	Description
325030	KRUSAN wall dispenser, 1000 ml
325031	KRUSAN wall dispenser, 1000 ml, short



KRUSAN Defiren Disinfection

The mixture is 1:50, which is corresponding to 1 dl KRUSAN Defiren mixed with 5 l water. The finished dissolution can be used for hands, arms, surfaces, instruments, inventory, etc. Can be used both indoors and outdoors.

Cat. No	Description
320425	KRUSAN DEFIREN disinfection and cleaning, 1000 ml, 6/pk





325141

325140

KRUSAN Chlorhexidine Conc

Effective in combating bacteria, fungi, spores and virus. Does not cross an intact skin barrier. Quick-acting and with long-term effect when used as a skin disinfectant. Chlorhexidine gluconate has continuous effects in the presence of organic substances, e.g. blood and saliva. For veterinary use in connection with obstetrics, dental work, wound treatment etc.

Cat. No	Description
325140	KRUSAN Chlorhexidine Conc., 5 %, with v chlorhexidine gluconate, 1 l
325141	KRUSAN Chlorhexidine Conc., 5 %, with v chlorhexidine gluconate, 2.5 l



325150



325151



325152

KRUSAN Chlorhexidine Alcohol

0.5 vol. % Chlorhexidine digluconate in 74 vol. % denat. alcohol. Ready mixed solution for disinfecting intact skin subsequent to washing. Application: to be used as a surgical hand disinfectant and as a disinfectant of the surgical site prior to operating. Store below 25°C.

Cat. No	Description
325150	KRUSAN Chlorhexidine alcohol 500 ml
325151	KRUSAN Chlorhexidine alcohol 1000 ml
325152	KRUSAN Chlorhexidine alcohol 2.5 l

**KRUSAN Brintoverilte**

En 3% vandig opløsning har en kraftig, men kortvarig desinficerende virkning. Opløsningen virker desuden vævssammentrækkende (adstringerende), idet protein udfældes i overfladen af det berørte væv, samt lugtdæmpende på grund af iltning af ildelugtende stoffer.

Cat. No	Description
325164	KRUSAN Brintoverilte 3% 5 ltr.





Obstetrics

To help our customers we have assembled the KRUUSE Obstetric Set including chain, ropes, handles, hook, calving belt as well as gloves and gel. Delivered in a sturdy and practical metal box.

You may also wish to have a look at some of the other products in the obstetrics range, such as our embryotome, uterine forceps and uterine pump.



KRUUSE Obstetric Set

Consists of:

- 180235: 2 x Nylon handles
- 180290: 1x Endless chain with 2 Harms eye hooks
- 180576: 1 x Pair of calving ropes 1 loop
- 180577: 1 x Pair of calving ropes 2 loops
- 180270: 2 x Eye hooks, blunt
- 180203: 1 x Calving belt with hook
- 180315: 1 x Obstetric hook Krey Schöttler
- 180426: 1 x Finger knife
- 180520: 2 x BOVIVET gel 500 ml
- 260700: 1 x KRUTEX Super Sensitive gloves
- 250170: 1 x Metal box with removable shelf and lid


Cat. No	Description
180600	KRUUSE obstetric set in metal box

Embryotome and Accessories

Thygesen model, nickel-plated with handle for fixation of obstetric chain.



Cat. No	Description
180003	Embryotome Thygesen, dismountable, 85 cm
180006	Embryotome Thygesen, 100 cm
180011	Embryotome Thygesen, 85 cm
180055	Jointed threader for Thygesen embryotome, 80 cm
180061	Jointed threader for Thygesen embryotome 1 m
180030	Tube end bush
180020	Spare head for embryotome
180040	Handle with groove for OB chain
180050	Clamping piece





Curved Snare Director

With screw clamp to hold saw wire.

Cat. No	Description
180120	Curved snare director

180146



Saw Wire Handle

The Danish saw wire handle and the Thygesen handle are made of stainless steel, while the squeezing type handle is chrome-plated and forged brass.

Cat. No	Description
180140	Danish saw wire handle
180146	Saw wire handle squeezing type
180145	Saw wire handle Thygesen



180140



180145



Saw Wire

Saw wire for embryotomy and dehorning.

Cat. No	Description
180150	BOVIVET saw wire 3.6 m, 12/pk
180165	BOVIVET saw wire 10 m
180160	BOVIVET saw wire 12 m
180170	JACO-VET saw wire 3.6 m
180171	JACO-VET saw wire 10 m

180160

180150

180171



180203

180202

BOVIVET Calving Belt

Calving belt with stainless hook. Simple and strong calving belt to be used together with obstetric chains or rope.

Cat. No	Description
180203	BOVIVET calving belt with hook
180202	BOVIVET Obstetric hook for calving belt



Obstetrical Chains

Made of stainless steel.

Cat. No	Description
180210	KRUUSE Obstetric chain 1 m, premium quality
180225	KRUUSE Obstetric chain 1.5 m, standard quality
180220	KRUUSE Obstetric chain 1.5 m, premium quality
180227	KRUUSE Obstetric chain with two handles, 1.3 m nickel-plated



KRUUSE Handles for Obstetric Chains

Cat. No	Description
180235	KRUUSE Nylon handle, 30 x 300 mm, with groove
180265	KRUUSE Handle with hook, 10 cm
180266	KRUUSE Handle with hook, 25 cm



Eye Hook

Made of stainless steel.

Cat. No	Description
180270	Harms eye hook blunt 8 cm
180275	Harms eye hook sharp 8 cm
180280	Harms obstetric stainless steel chain with single hook



Calfpullers

Cat. No	Description
180565	HK calf puller E11 with 2 angle pieces
180570	Calf puller
180575	Rope for calf puller model E11 2/pk
180571	Calving rope/helper (Vink)



KRUUSE Obstetric Hooks

Cat. No	Description
180310	Krey-Shottler obstetric hook, sharp
180315	Krey-Schottler obstetric hook, standard



KRUUSE Calving Rope

Cat. No	Description
180576	Calving rope with 1 loop pair, Ø=5 mm, length 150 cm
180577	Calving rope with 2 loops pair, Ø=5 mm, length 150 cm




BOVIVET Vaginal Suture Tape


Cat. No	Description
152130	BOVIVET vaginal tape 5 m roll
152135	BOVIVET vaginal tape 0.8 m sterile, 12/pk

Vulva Suture Needles

Complete with hose as well as large and small connector.



Cat. No	Description
152110	Gerlach suture needle, stainless steel, 19 cm
152120	Bühner suture needle, stainless steel, 30 cm



BOVIVET Prolapse Pins

With plastic or wooden finger grip, adjustable or fixed.

Cat. No	Description
152145	BOVIVET prolapse pins 55-75 mm, 12/pk
152149	BOVIVET prolapse pins with wooden balls 40 mm, 12/pk
152150	BOVIVET prolapse pins with wooden balls 50 mm, 12/pk
152151	BOVIVET prolapse pins with wooden balls 60 mm, 12/pk
152152	BOVIVET prolapse pins with wooden balls 65 mm, 12/pk
152153	BOVIVET prolapse pins with wooden balls 70 mm, 12/pk
152160	Plastic splints for prolapse pins, 12 pairs
152170	Plastic splints for prolapse pins with thread, 10 x 1 pair



BOVIVET Teat Wipe

For disinfecting the teat before treatment.

Cat. No	Description
170100	BOVIVET teat wipe 100/pk



BOVIVET Silicone Teat Insert

The Bovivet silicone teat insert is a gentle teat insert designed with a handle for easy and hygienic insertion. The conical shape of the teat insert prevents the teat insert from falling out.

Cat. No	Description
170230	BOVIVET silicone teat insert, sterile, 20/pk



BOVIVET Teat Dilator

Diam. 5.5 mm, length 35 mm, individually packed in plastic tube with vaseline, red obturator to be removed when milking.
For cows with damaged teats, the drain can be inserted for an extended period.

Cat. No	Description
170270	BOVIVET teat dilators, 5 x 10/pk
170275	BOVIVET teat dilators, 25 x 10/pk
170280	BOVIVET teat dilators thin, 5 x 10/pk



BOVIVET Wax Teat Plugs

Smooth and flexible wax teat plug with head for easy and sterile application. This teat plug will slowly dissolve in the teat canal. Recommended by leading specialists.

Cat. No	Description
170283	BOVIVET wax teat plugs with head, sterile, 20/pk



BOVIVET Teat Plugs Long

For single use, extra long, very flexible. Packaged in a 5 cm diameter box with vaseline. Universal connector. Red pin can be removed before milking.

Cat. No	Description
170290	BOVIVET teat plugs long, 7 cm, 10/pk



KRUUSE Colombus Teat Plugs in Box

Single use. Box with vaseline. Each plug placed in its separate chamber.

Cat. No	Description
170300	KRUUSE Colombus teat plugs in plastic box, 20/pk



BOVIVET Teat Dilators

Single use. Bengen model, diameter 2.5 mm, length 4.5 cm, with vaseline. Use to open blocked teat channel.

Cat. No	Description
170315	BOVIVET teat dilators, 6 x 20/pk



BOVIVET Teat bandage

Unique polyester bandage. Protects wound against dirt, easy to mold and apply, cohesive. Can be applied to wet teats. Allergy tested.

Cat. No	Description
160954	BOVIVET Teat bandage, 6 cm x 4.5 m



BOVIVET Milking Tube

Yellow, 3.3 mm. Sterile milking tube in blister pack. Luer lock for easy administration of e.g. antibiotics.

Cat. No	Description
170193	BOVIVET milking tube luer lock, yellow, 3.3 mm, 20/pk, sterile



BOVIVET Syringe Mount Strips

Used for administration of medication from a syringe into the teat. Luer fitting.

Cat. No	Description
170260	BOVIVET syringe mount strips, luer, 100/pk, sterile



Milk Catheter

Complete with stiletto, with olive hub, nickel-plated.



Cat. No	Description
170165	Milk catheter 2.0 mm x 9 cm
170170	Milk catheter 2.5 mm x 9 cm
170180	Milk catheter 3.0 mm x 9 cm
170190	Milk catheter 4.0 mm x 9 cm

KRUUSE Teat Drainage Tubes

With stylet and teat obturator.



Cat. No	Description
170250	KRUUSE teat drainage tubes, 7 cm, 25/pk

KRUUSE Teat Dipper

For disinfection of the teats before milking. Consists of two chambers. Pressure on the lower chamber forces the liquid into the upper cup, and the teat is easily washed.



Cat. No	Description
170339	KRUUSE teat dipper - non return valve

BOVIVET Milk Sucking Preventer

Aluminium, with spikes.



Cat. No	Description
170380	BOVIVET sucking preventer spiked

KRUUSE Plastic Sucking Preventer

Nylon, made of flexible plastic. Easy to apply.



170385

170386

Cat. No	Description
170385	KRUUSE Plastic sucking preventer, S
170386	KRUUSE Plastic sucking preventer, M



Cattle dehorning is a normal part of management practice to prevent injury to animals and handlers.

It is generally accepted that the younger the calf is at the time of dehorning, the less invasive and painful the procedure. KRUUSE offers gas and electric dehorning devices. Higher working temperatures can be obtained using the gas dehorning devices. We have dehorning devices for calves as well as for adult animals.



BOVIVET Electric Dehorning

Very solid and reliable dehorning device including heating element enabling the brass tip to reach a working temperature of 700° C. Supplied with 2 m cord.

Cat. No	Description
240000	BOVIVET electric dehorner 220 V complete with holder and 2 dehorning tips
240001	BOVIVET electric dehorner, 230 V, complete with holder and 2 dehorning tips, CE approved
240003	Dehorning tip, sharp, for BOVIVET electric dehorner
240007	Dehorning tip, flat for BOVIVET electric dehorner



Gas Dehorning

The dehorner MR III's working head has increased strength with the addition of a 20 mm stainless steel tube flame housing. Resulting in the dehorner III's front end being much stronger and robust.

Cat. No	Description
240260	Dehorner MR III with electronic ignition
240261	Dehorning spare tip 15 mm with catalyzer for 240260
240262	Dehorning spare tip 18.5 mm with catalyzer for 240260
210096	KRUUSE Lighter gas 200 ml for Dehorner



BOVIVET Anti Kick Hobble

With variable release catches, to keep hind legs in place after calving or milk fever.

Cat. No	Description
210058	BOVIVET nylon anti kick hobble



Styrko Kick Stopper

Extra strong, adjustable with spring lock, prevents kicking at milking and udder treatment.

Cat. No	Description
210060	Styrko cow kick stopper
210080	Replacement springs - D 9 mm



Bull Holder

Cat. No	Description
210315	Bull holder, Harms model standard quality, with rolling adjuster, 19 cm
210316	Bull holder, Harms model standard quality, with slide adjuster, 19 cm



Nose Clamp

With spring. Easy to attach and remove.


Cat. No	Description
210380	Nose clamp for bull 11 cm



Nose Clamp

Bull holder, 100 cm, American model.

Cat. No	Description
210325	Nose clamp for cattle with nylon cord, standard quality



Bull Rings

With security screw. Standard quality. Stainless steel.

Cat. No	Description
210515	Bull ring 55 mm, stainless steel, 6/pk
210516	Bull ring 59 mm, 6/pk
210518	Bull ring 55 mm, brass, 6/pk
210540	Bull ring 65 mm
210545	Bull ring 65 mm, chrome plated



Bull Ring Applicator

With turning handle for inserting 55 mm bull rings. Length 340 mm.

Cat. No	Description
210450	Bull ring applicator for 55 mm rings





KRUUSE PASTE

Recovin – a series of products for oral administration of specific nutrients for primarily dairy cows.

Within the Recovin range you will also find pastes for other species such as sheep and goats.

Paste administration is easy and safe with the Recovin dosing gun and only takes 30 seconds. The full range has been developed and tested by Danish veterinarians and all products are made in Denmark.



ReCovin Ketosis Paste

A high energy product containing propylene glycol and other relevant essential nutrients for ruminants in negative energy balance.

Cat. No	Description
310760	ReCovin ketosis paste 4 x 250 g



ReCovin Rumen Acidosis Paste

Acidosis paste contains magnesium oxide with antacid effect and other useful nutrients for ruminants with rumen acidosis.

Cat. No	Description
310764	ReCovin rumen acidosis paste 4 x 280 g



ReCovin Calcium Paste

Recovin calcium paste is a high-tech calcium chloride product with added magnesium. Recovin Calcium Paste is especially useful as a calcium supplement for cows in the period around calving where they are at risk of milk fever.

Cat. No	Description
310766	ReCovin calcium paste 4 x 295 g



ReCovin Rumen Starter

ReCovin rumen starter paste is a unique paste, with nutrients that act as a substrate for the build up of the microflora in the rumen after periods of ruminal disturbances and decreased appetite. ReCovin rumen starter contains vitamins, micro minerals (trace elements), amino acids, energy and appetizing agents.

Cat. No	Description
310773	ReCovin rumen starter 4 x 255 g



ReCovin Phosphate Paste

ReCovin phosphate paste is used as a source of phosphate for cows with downer cow syndrome (cows with low serum organic phosphorus that do not respond to intravenous calcium therapy). Phosphate paste is phosphate mixed with the energy provider propylene glycol.

Cat. No	Description
310880	ReCovin phosphate paste 4 x 350 g



ReCovin ACT Paste

Active coal in form of pasta, contains tannin. For ruminants and horses.

Cat. No	Description
310767	ReCovin ACT paste 4 x 250g



ReCovin Dosing Gun

ReCovin paste dosing gun. Very durable. This model does not have variable dosing facility. One dose equals 10 ml pasta.

Cat. No	Description
311010	ReCovin dosing gun, plastic and metal (black)



Artificial Insemination - Bovine

The KRUUSE bovine AI range includes commonly used instruments as well as catheters, sheaths and gel.

To make ordering easy, we have assembled the KRUUSE Artificial Insemination Kit which includes many of the items from this chapter as well as disposable gloves, apron, scissors and a practical and sturdy stainless steel box.



KRUUSE Artificial Insemination Kit

- Set containing:
- 1 x Inseminator universal
 - 1 x Pack of sheaths universal, 50/pk
 - 1 x Aluminium tube for sheaths
 - 1 x Straw tweezer
 - 1 x Pair of scissors 14 cm
 - 1 x Vaginal speculum with light
 - 1 x Thermo container 0,5 l
 - 4 x Bottles of gel sterile 500 ml
 - 1 x KRUTEX Super Sensitive disposable gloves, 100/pk
 - 1 x KRUTEX obstetric apron 120 cm
 - 1 x Metal box stainless steel

Cat. No	Description
340680	Artificial insemination kit for cattle

Artificial Vagina Bull

Complete with rubber liner, collecting cone and collecting tube.
The vagina can be inflated with air or water.



Cat. No	Description
340284	Artificial vagina for bull incl. Rubber cone and collecting tube
340330	Rubber liner f/bull and boar vagina Ø80 mm x 50 cm, smooth
340350	Collecting tube 15 ml
340600	Spare cone for artificial vagina





Vaginal Speculum Polansky

The speculum has 3 blades. Has butterfly screws with 2 blades on top for wide exposure.

Nickel plated

Cat. No	Description
290990	Vaginal speculum Polansky 27 cm
290995	Vaginal speculum Polansky 35 cm



Vaginal Speculum

Vaginal speculum for horses and cows. Nickel plated.

Cat. No	Description
290985	Vaginal speculum, 27 cm metal duck bill type
290980	Vaginal speculum, 35 cm metal duck bill type



KRUUSE Gel, Sterile

Excellent lubrication for AI catheters and instruments. Non-spermicidal. No preservatives.

Cat. No	Description
340623	KRUUSE sterile gel 10 ml
340625	KRUUSE Sterile gel 250 ml

**KRUUSE Catheter Gel**

Efficient lubricating gel to be used e.g. in connection with the introduction of AI catheters. Handy bottle.


- Non-spermicidal
- No preservatives
- Contents: Vaseline and paraffin oil

Cat. No	Description
340626	KRUUSE catheter gel, hanging tube, 250 ml

**Bovine Uterine Metal Catheter**

Bovine uterine metal catheter 43 cm, luer lock

Cat. No	Description
230710	Bovine uterine metal 43 cm, catheter luer lock

**Inseminator Universal**

For use with mini and medium straws.

Cat. No	Description
340665	Inseminator universal for both 0.25 ml and 0.50 ml straws

**Aluminium tube for sheath**

Sheat holder in aluminium

Cat. No	Description
340660	Sheat holder in aluminium

Sheaths for Universal Inseminator

For insemination gun. Fits mini and medium straws.



Cat. No	Description
340655	Sheaths, split, universal for insemination gun, 50/pk
340656	Sheaths, non-split, universal for insemination gun, 50/pk

Straw Tweezers

The stainless steel straw tweezers secure easy handling of mini and medium straws.



Cat. No	Description
340315	Straw tweezers, 15 cm, stainless steel

Straw Cutter

Designed to quickly cut off the sealed end of straws. After insertion of straw a push on the button cuts the straw at the proper length.



Cat. No	Description
340320	Straw cutter

KRUUSE Insemination Catheter




Cat. No	Description
230830	KRUUSE Insemination Catheter 53.5 cm, sterile, 25/pk
230831	KRUUSE Insemination Catheter 60 cm, sterile, 25/pk
230832	KRUUSE Insemination Catheter 70 cm, sterile, 25/pk
230840	KRUUSE Insemination Catheter 53.5 cm, non-sterile, 25/pk
230841	KRUUSE Insemination Catheter 60 cm, non-sterile, 25/pk
230842	KRUUSE Insemination Catheter 70 cm, non-sterile, 25/pk



Kamar Heat Detector

Helps identifying standing heats. Catches cows at night when majority of standing heats occur. Lowers per heat insemination expenses.

Cat. No	Description
340645	Kamar heat detector 25/pk



Nitrogen Container Forceps

For removal of goblets from nitrogen container, stainless

Cat. No	Description
340290	Nitrogen container forceps 70 cm
340620	Nitrogen container forceps 25 cm



Thawing Unit for Straws

The unit is electronically controlled for thawing frozen semen. It automatically heats and maintains an accurate water temperature at 35,6°C. Can be used for mini and medium type of straws. Supplied with car lighter plug, 12 V.

Cat. No	Description
340174	Thawing unit 12-220V for straws 35-43 degr. C

340127



340131

Semen Container GT

Nitrogen container for storage and transport of semen.

Cat. No	Description
340127	Container for semen GT2, 2 l Max. capacity of 0.25 ml straws - 330, max. capacity of 0.5 ml straws - 50. Static holding time 25 days, dynamic holding time 15 days.
340131	Container for semen GT 9, 9.3 l Max. capacity of 0.25 ml straws - 1560, max. capacity of 0.5 ml straws - 720. Static holding time 84 days, dynamic holding time 50 days.
340133	Container for semen GT11, 12.2 l Max. capacity of 0.25 ml straws - 1560, max. capacity of 0.5 ml straws - 720. Static holding time 130 days, dynamic holding time 80 days.
340136	Container for semen GT21, 21.5 l Max. capacity of 0.25 ml straws - 1560, max. capacity of 0.5 ml straws - 720. Static holding time 225 days, dynamic holding time 140 days.
340141	Container for semen GT35, 33.6 l Max. capacity of 0.25 ml straws - 7380, max. capacity of 0.5 ml straws - 3285. Static holding time 90 days, dynamic holding time 56 days.
340150	Container for semen GT40, 40 l Max. capacity of 0.25 ml straws - 16400, max. capacity of 0.5 ml straws - 7300. Static holding time 140 days, dynamic holding time 88 days.



340147



340149

Nitrogen Container TR

Container for storage and transport of nitrogen.

Cat. No	Description
340147	Container for transport TR7, 7,2 l Static holding time 36 days.
340148	Container for nitrogen TR11, 12.2 l Static holding time 67 days.
340149	Container for nitrogen TR26, 26.0 l Static holding time 130 days.

Semen Container Voyageur

Nitrogen container especially designed for difficult transports.

The containers are equipped with a material which absorbs the liquid nitrogen and therefore the containers cannot leak.



340126



340137

Cat. No	Description
340126	Nitrogen container Voyageur 2, 1.75 l Max. capacity of 0.25 ml straws - 220, max. capacity of 0.5 ml straws - 100. Static holding time 13 days, dynamic holding time 8 days.
340137	Nitrogen container Voyageur 5, 6.5 l Max. capacity of 0.25 ml straws - 1040, max. capacity of 0.5 ml straws - 400. Static holding time 37 days, dynamic holding time 23 days.



Semen filters

Cat. No	Description
340171	KRUUSE Gauze pads 30x30 cm, 360/pk
340173	Ultra-flow boar semen filter Ø22.9 cm, 200/pk



Water Bath, Semen Diluent

Plastic water bath for reheating of semen diluent, 12 l, 850 W, adjustable temperature, digital display.

Cat. No	Description
340120	Water bath for reheating semen diluent
340125	Water bath 28 L

Laboratory Thermometer

Cat. No	Description
291125	Laboratory thermometer-10/+110°C. Length 26 cm.





Mobile Storage Box Boar Semen

Ideal for the car or when high mobility on the farm is needed.

Integrated heating and cooling system.

Temperature on digital control, preset at 17°C.

Power supply 220 V, 12 V DC, 24 V DC.

Power consumption as low as 64 W Optimal temperature control when storing boar semen is important. The KRUUSE range of storage equipment for boar semen ensures that valuable semen is stored under optimal conditions, ie at a temperature of 16-18°C. High quality parts make up this exceptional storage equipment, and all devices have undergone 3 days of thorough testing at the factory prior to dispatch in order to ensure high performance and uncompromising quality.

Cat. No	Description
340122	Climate control box 35 l for transport of boar semen +18C, 220V/12V. WxDxH: 390x560x500 mm



Storage Cabinet Boar Semen

Integrated heating and cooling system.

Ventilation for optimal temperature control and constant temperature for the entire unit.

Power supply 220V.

Power consumption as low as 100 W.

Adjustable temperature on digital control, preset at 17°C Optimal temperature control when storing boar semen is important. The KRUUSE range of storage equipment for boar semen ensures that valuable semen is stored under optimal conditions, ie at a temperature of 16-18°C. High quality parts make up this exceptional storage equipment, and all devices have undergone 3 days of thorough testing at the factory prior to dispatch in order to ensure high performance and uncompromising quality.

Cat. No	Description
340112	Climate control cabinet for boar semen, 130 l



KRUUSE BTS Plus boar semen diluent

KRUUSE Boar semen extender, mixed with distilled water

Ensures 3-5 days semen conservation

Easy and quick preparation of diluent

Contains gentamycin and marbofloxacin

Cat. No	Description
340020	KRUUSE BTS Plus boar semen diluent, 1 l pouch
340021	KRUUSE BTS Plus boar semen diluent, 5 l pouch



KRUUSE Durasperm Plus boar semen diluent

KRUUSE Boar semen extender, mixed with 1 liter of distilled water
Ensures 5-7 days semen conservation
Easy and quick preparation of diluent
Contains gentamycin and marbofloxacin

Cat. No	Description
340023	KRUUSE Durasperm Plus boar semen diluent, 1 l pouch
340024	KRUUSE Durasperm Plus boar semen diluent, 5 l pouch



Demineralized water

Demineralized sterile-filtered water, 5 l. Shelf life of 6 months from date of production in original undamaged packaging.

Cat. No	Description
340169	Demineralized sterile-filtered water 5 l



Semen Bags

Cat. No	Description
340519	Semen bags, 3 l, 150/pk
340523	Semen bags, 5 l, 50/pk

Semen Bottles

Cat. No	Description
340630	Semen bottle, 100 ml without lid, 100/pk, (suitable for lids: cat. No 340631-340635)
340631	Lid for semen bottle (cat. No 340630), blue, 100/pk
340632	Lid for semen bottle (cat. No 340630), yellow, 100/pk
340633	Lid for semen bottle (cat. No 340630), green, 100/pk
340634	Lid for semen bottle (cat. No 340630), green, 100/pk
340635	Lid for semen bottle (cat. No 340630), green, 100/pk



KRUUSE Insemination Catheters

KRUUSE A.I. catheters have rounded edges preventing penetration of the semen bag for hygienic and easy introduction in the sow.

Total length: 550 mm
Outer diameter: 7.0 mm
Inner diameter: 4.9 mm
Length of foam tip: 32 mm

Cat. No	Description
340210	KRUUSE foam-tip insemination catheter, rounded edge, 50 x 10/pk
340208	KRUUSE insemination catheter w/ yellow foam tip, rounded edge, single packed in plastic sheath



KRUUSE Lubricated Catheters

Single packed, securing optimal hygiene conditions and improvement of conception rate. The advantages of clean insemination are well known to pig breeders. The avoidance of bacterial contamination of the uterus during insemination improves the performance in terms of conception rate and litter size. The protective and hygienic foil is used as a protection sheath during insertion of the catheter. The foam-tip catheter is lubricated with the well-known KRUUSE catheter lubricating gel, tested by the Danish Meat Association and approved to be non-spermicidal.

Cat. No	Description
340202	KRUUSE lubricated catheter, single packed, 100/pk

KRUUSE Insemination Catheter with Handle / Stopper to prevent Back Flow of the Semen

Specifications: Total length 535 mm and length of foam tip 32 mm.



Cat. No	Description
340214	KRUUSE foam-tip insemination catheter, with handle, 100 x 5/pk, 20/pk



KRUUSE Insemination Catheter for Gilts

KRUUSE A.I. catheters have rounded edges preventing penetration of the semen bag for hygienic and easy introduction in the sow.

KRUUSE insemination catheter with small foam tips suitable for gilts.

Total length: 550 mm

Outer diameter: 7.0 mm

Inner diameter: 4.9 mm

Length of foam tip: 34 mm

Cat. No	Description
340220	KRUUSE foam-tip gilt insemination catheter, rounded edge, 100 x 5/pk, 20/pk




KRUUSE Spiral Insemination Catheters

KRUUSE disposable spiral catheters are soft and gentle. They are delivered with a tube connector with closure mechanism preventing backflow.

Made of hygienic material.

Cat. No	Description
340198	KRUUSE spiral insemination catheter, with handle, 4 x 1/pk, 125/pk



KRUUSE intra-uterine insemination catheter

Cat. No	Description
340190	KRUUSE intra-uterine insemination catheter



KRUUSE Liquid Catheter Gel

Efficient lubricating gel to be used e.g. in connection with the introduction of AI catheters. Handy bottle. Non-spermicidal. No preservatives. Contents: Vaseline, paraffin oil.

Cat. No	Description
340626	KRUUSE catheter gel, hanging tube, 250 ml

**KRUUSE Sterile Gel**

Excellent lubrication for AI catheters and instruments. Non-spermicidal. No preservatives.

Cat. No	Description
340625	KRUUSE Sterile Gel 250 ml, 18/pk

**PORCIVET Boar Spray**

PORCIVET boar spray is used for providing standing oestrus during insemination of sows and gilts. CO₂ is used as propellant. From a distance of approx. 60 cm spray the Porcivet boar spray directly onto the snout of the sow/gilt. Spray for approx. one second and carry out the standing test. If the sow willingly accepts the pressure, she is ready for insemination.

Cat. No	Description
241129	PORCIVET Boar Spray 110 ml, 12/pk
241124	PORCIVET Boar Spray 250 ml, 12/pk

Ultrasound Scanner

The hand held WED-2000A is distinguishable by the light weight, easy hand carrying, accurate imaging and simple interface. The selectable dual frequency probe offers a more widely range of applications.

The WED-2000A offers a cost effective solution to vets, A. I. technicians and pig farmers who need to carry out quick and easy diagnosis for pregnancy detection and sexing. Technical specifications: Monitor: 5.0 inch TFT LCD Weight: 750 g

Probes: 3.5 MHz/5.0MHz Mechanical Sector Selectable

Battery: Rechargeable and internal, back-up for 4hrs

Display Mode: B,B+B,B+M,M

Image magnification in four levels:

1.0, 1.2, 1.5, 2.0 times

Gray Scale: 256 levels

Text Annotation: Patient, ID, sex, age

Software: For human or veterinary, selectable

Measurement:

Distance, length, area/circumference, time/heart rate/EF rate for heart

Stander Accessories: AC Adapter, internal battery, gel, mouse

Optional: 5.0 MHz sector probe, car charger

Extra battery with charging cable

Standard configuration: main unit + 3.5 MHz waterproof probe

Options:

1). 3.5 MHz waterproof probe

2). Extral rechargeable internal battery with the cable

3). Car Charger



Cat. No	Description
241140	WED-2000AV Ultrasound scanner for pigs



MSU 1 Ultrasound Scanner
Full-digital Mechanical Sector Ultrasound Scanner

Function feature

- Display mode: B, 2B, 4B, B/M, M.
- Image storage: 100 frames permanent storage, without loss of image when powered off.
- Image Processing: IP (4 stages of Gamma Correction), Frame correlation processing (4 stages), Left/Right, Up/Down reverse.
- 5 kinds of magnifications.
- Interface language: English/Chinese
- Internal pseudo color display (8 kinds).
- Body mark: 27 kinds including Swine, Ovine, Canine, Feline, Goat, Llama ect.
- High-speed USB socket, support large-capacity flash disk and mouse.
- Puncture guide and gravel positioning function.

Measurement function

- Common measurement: distance, circumference, area (tracing method, elliptic method), volume, ect.
- OB measurement: Bovine, Equine, Ovine, Canine, Feline, Goat, Llama, gestational weeks and expected date.
- Cardio measurement: depth, slope, heart rate, cycle.

Report function

- OB report

Cat. No	Description
241201	Full-digital Mechanical Sector Ultrasound Scanner



Image quality



Gravidity of pig which is 25 days



Gravidity of pig which is 86 days



Gravidity of dog which is 35 days



Dog's bladder

MSU 2 Ultrasound Scanner

Full Digital B Mode Veterinary Ultrasound Scanner

Function feature

- Display mode: B, 2B, 4B, B/M, M, total 5 kinds
- 100 frames permanent storage
- Local zooming
- Eight internal pseudo colors (incl. B/W)
- OB measurement: EDD and GA for swine, ovine, and canine ect.
- OB report
- Flash card storage
- Dual power supply, equipped with high capacity lithium battery, driving the unit more than 6 hours.

Pregnancy diagnosis

- MSU 2 can check whether the pig is pregnant 18 days later mating or not.
Comparing with traditional method, you can get non-pregnant pigs about 12 days earlier.

Diagnosis towards fetus

- Measure heart size and watch for heartbeat.

Diagnosis of sow uterine function

- Diagnose sow reproductive disturbance diseases: Pyometra, mucometra, etc.

Standard configuration:

- Main unit: 1 PC
- 3.5 MHz mechanical sector probe: 1 PC
- Lithium battery: 1 PC
- Base Changer: 1 PC
- Sunshine cover: 1 PC
- Silicone Protective Cover: 1 PC
- Plastic sealed box: 1 PC

Optional configuration:

- Lithium battery
- Black leather sabbiebag
- Goggles

Cat. No	Description
241202	Full Digital B Mode Veterinary Ultrasound Scanner

Image quality



Non-pregnant uterine horn



Gravidity of pig which is 28 days



Mucometra and pyometra



Dog's bladder



Echophonic Scanner Gel

The ultimate gel for advanced ultrasound examinations. This gel is water soluble and hypoallergenic. Tissue and transducer friendly. After use remove the gel with tissue or rinse with water. Optimum viscosity and minimum air content ensures the best possible contact for optimum performance of your ultrasonic probes. 5 l can be delivered with an empty 250 ml application bottle.

Cat. No	Description
290690	Echophonic ultrasound transmission gel, 250 ml, 12/pk
290692	Echophonic ultrasound transmission gel, 5 l



PORCIVET Shoulder Protector

Developed to prevent shoulder wounds on sows. The shoulder cap is designed to relieve pressure in the area around the shoulder blade. This effectively reduces the occurrence of shoulder wounds and subsequent treatments or destruction of the animal. The protector uses a unique adhesive that can be adhered directly to the affected area of the sow.


Cat. No	Description
241013	PORCIVET shoulder protector, 10/pk



KRUUSE Coaxer

Handy pocket model, solid, length 15 cm, uses 3 x 1.5 V batteries.

Cat. No	Description
210045	KRUUSE coaxer without battery



Hot-Shot POWER-MITE

HOT SHOT coaxer, Power-Mite Compact, is a very strong model, measures 15.6 x 5 x 3.2 cm, is designed for close-in work, ideal for veterinary surgeons to carry in their bag.

Cat. No	Description
210040	Hot-Shot Power- mite coaxer without battery




Hot-Shot Dura Prod

New and improved design with 81 cm flexible shaft, yellow.

Output 9.000 V at 15 – 36 Watts.

Features include: sealed circuitry, patented technology to reduce feedback shock.


Cat. No	Description
210026	Hot-Shot flexible coaxer, with battery, yellow
210029	Flexible extension handle for Hot-Shot coaxer 32"/81 cm



Hot-Shot Coaxer HS 2000

With 81 cm flexible shaft, providing the highest output from the least amount of battery draw. Offers the highest performance in all weather conditions, especially on wet animals. Sealed motor technology, cannot be damaged by dirt or moisture. Maximum spark and battery life, easy access to battery and motor compartment, uses 6 x 1.5 V batteries.

Cat. No	Description
210030	Hot-Shot flexible coaxer with battery, green
210031	Flexible extension handle for Hot-Shot coaxer, 42"/107 cm
210036	Motor for green Hot-Shot coaxer



Hot-shot Coaxer Sabre-Six

With 81 cm flexible shaft and effective medium power output. Features the sealed motor technology which means that the electronics inside the unit are protected and will not be damaged by dirt, manure, and moisture. Works even on wet animals. Easy access to battery and motor compartment, high performance at a reasonable price, uses 6 x 1.5 V batteries.

Cat. No	Description
210032	Hot-Shot flexible coaxer, with battery, red, 6/pk
210033	Flexible extension handle for Hot-Shot coaxer, 22"/56 cm
210037	Motor for red Hot-Shot coaxer





KRUUSE Heating Lamp

Long life light bulb of 5,000 hours. Shatter resistant of hard glass. Splash proof with incandescent reflection. With threaded laminated metal base, which means you never lose base. Aluminium insulator 220-240 V.

Cat. No	Description
210140	KRUUSE heating lamp, 100W, clear, 2/pk
210141	KRUUSE heating lamp, 150W, clear, 2/pk
210143	KRUUSE infrared heating lamp, 100W, red, 2/pk
210144	KRUUSE infrared heating lamp, 150W, red, 2/pk
210145	KRUUSE infrared heating lamp, 250W, red, 2/pk
210152	KRUUSE PAR-lamp heating lamp, 100 W, red, 1/pk
210153	KRUUSE PAR-lamp heating lamp, 100 W, clear, 1/pk



Sorting Plates

Very robust plate for sorting livestock.

Cat. No	Description
210006	Pushing plate, red, 46 x 76 cm
210007	Pushing plate, red, 96 x 76 cm
210008	Pushing plate, red, 120 x 76 cm



Sorting Paddle

Another humane alternative for livestock sorting and handling. The brightly coloured paddle head (31 cm long and 16 cm wide) acts as a visual guide and the small ball bearing pieces inside make a sound to which animals respond. High quality. Vinyl covering for safety and easy cleaning, quality golf-style handle for a sure grip and easy manoeuvre, an excellent alternative when you do not want to use a coaxer.

Cat. No	Description
210011	Sorting paddle 107 cm, red
210012	Sorting paddle 122 cm, blue



KRUUSE Sorting Paddle


Another humane alternative for livestock sorting and handling. The brightly coloured paddle head (31 cm long and 16 cm wide) acts as a visual guide and the small ball bearing pieces inside make a sound to which animals respond.

Cat. No	Description
210014	KRUUSE Sorting paddle, non-electric, 107 cm



Pig Holder

Cat. No	Description
210340	Pig holder, steel cable 25 cm
210344	Pig holder, Trops for piglets less than 30 kg
210343	Pig holder, Trops with steel wire
210346	Pig holder with handle, rigid model



Antonius Farrowing Snare

Cat. No	Description
180405	Antonius farrowing snare, length 560 mm
180407	Thread stainless for Antonius farrowing snare
180420	KRUUSE Obstetric hook for pigs blunt

KRUUSE Tooth Grinder

Professional tooth grinding machine especially developed for the grinding of piglets' teeth. Works approx. 8 hours per charge.



Cat. No	Description
240800	Tooth grinder International model incl. charger, battery, grindstone and plastic cover
240805	Grindstone for tooth grinder
240807	Charger for international tooth grinder
240808	Battery for International tooth grinder
240809	Plastic cover for the tooth grinder
240806	Motor for International tooth grinder





KRUUSE Pig Tail Docker

Docks and cauterizes in one action reducing risk of infection. Tail docking should be carried out shortly after birth (max. 2-4 days). The working temperature of the knife is 400°C. Works with standard lighter gas. All metal parts are stainless steel. Economical - docks up to 120 piglets per refill. The complete set consists of one tail docker, one bottle of 200 ml lighter gas and instructions for use in a durable box.

Cat. No	Description
210095	KRUUSE pig tail docker complete
210094	KRUUSE pig tail docker without gas
210102	Empty gas cartridge without knife for tail docker
210100	Spare knife for tail docker



Super Pro Pig Tail Dockers

For the professional pig breeder who requires high capacity and reliable operation. The working temperature of the knife is 500-525°C. Requires only 2 minutes for pre-heating. Piezo-ignition and works with standard lighter gas.

All metal parts are stainless steel. Economical - docks up to 150 piglets per refill. The complete set consists of one tail docker, one bottle of 200 ml lighter gas and instructions for use in a durable box.

Cat. No	Description
210088	Super Pro tail docker with gas
210089	Super Pro tail docker without gas
210103	Gas burner with knife for Super Pro tail docker
210101	Spare knife for pig tail docker SuperPro



Lighter Gas for Tail Docker

Cat. No	Description
210096	Lighter gas 200 ml for tail docker

240570



Tooth Cutting Forceps

Both models are with spring.

Cat. No	Description
240566	Tooth cutting forceps with short jaws 13 cm, stainless steel
240570	Tooth cutting forceps with long jaws 15 cm, stainless steel

240550



240560



240551



Ear Notcher

Very sharp cutting that gives a clean cut in the ear.

Cat. No	Description
240550	Ear notcher for calves, 16 cm, stainless steel, U-shape, 5 mm large
240560	Ear notcher for pigs, 16 cm, stainless steel, V-shape
240551	Ear notcher for piglets, U-shape, 4.5 mm large

KRUUSE Castration Forceps Pigs

Joins forceps and scalpel in one tool. Standard scalpel blades apply, no. 24. Secures lower rate of ruptured piglets. Ergonomic and safe design. Manufactured from stainless steel.

Cat. No	Description
240574	KRUUSE combi castration forceps
141182	KRUUSE sterile scalpel blade no 24, stainless steel, 100/pk, 20/pk



Omniplast

A hypoallergenic tape, adhering to the skin without risking allergic reactions. Easy to tear off, and can be used to fix medical instruments and dressings on necrotic tissue.

Cat. No	Description
169006	Omniplast 2.5 cm x 5 m, 12/pk
169007	Omniplast 5 cm x 5 m, 6/pk





Tally Counter

An easy way to count your animals. 4 digits, counts 0-9999, chrome-plated.

Cat. No	Description
247135	Tally counter



Umbilical Cord Clamps

Disposable, 65 mm.

Cat. No	Description
142131	Umbilical cord clamps, 100/pk



Obstetric Gloves Pigs

This glove is for assisting in the birth of piglets, particularly useful when a piglet is stuck. Obviates the need for using instruments.

Cat. No	Description
210109	Obstetric gloves for pigs, 5/pk



O-Bite and Hartshorn Spray

O-bite pig perfume is very effective against tale bits and pigfights. O-bite can also be used on chicken and poultry. PORCIVET Hartshorn oil discourages fighting and cannibalism among pigs and poultry.

Cat. No	Description
161072	O-bite spray, 400 ml, 12/pk
161073	PORCIVET hartshorn oil, 400 ml, 6/pk



Feeding Tubes

To be used for tube feeding of piglets. Metric tubes have centimetre graduation starting 5 cm from the tip of the tube. Easy to determine the accurate position of the tube. Frosted surface makes lubricant superfluous. Connector fits luer syringes and syringe pumps. Connector with plug provides hygienic protection against reflux and contamination. Radiopaque line enables checking the position of the tube. The tube can be used for up to 4 weeks without losing its softness.

Cat. No	Description
230800	Feeding tube for pigs CH 8 sterile, 10/pk, 100/pk
230801	Feeding tube for piglets CH10, sterile, 10/pk, 100/pk



KRUUSE Ovi-mark Marking Sprays

Special marking spray for sheep. Colours remain visible 5-6 months. Fully scorable. Approved by The British Wool Marketing Board.

Cat. No	Description
240460	KRUUSE Ovi-mark, red, 400 ml, 12/pk
240461	KRUUSE Ovi-mark, blue, 400 ml, 12/pk
240462	KRUUSE Ovi-mark, green, 400 ml, 12/pk



Sheep Restrainer

Unique plastic restrainer to aid the shepherd in his everyday tasks. Easy to fit, causes the ewe minimum discomfort. Ideal for lambing, transport and during treatment.

Cat. No	Description
247090	Sheep restrainer plastic large



Lambing Instrument

Useful tool for assisting lambing. Snare can be put over head or leg, pulled tightly, and the lamb is then pulled out gently. Plastic, length 37 cm.

Cat. No	Description
247080	Lambing instrument plastic



Uterine Support

This plastic support is used to prevent prolapse around the time of lambing. Once inserted into the vagina it can be tied to the fleece. Sterilize and use a lubricant before insertion.


Cat. No	Description
247000	Uterine support, plastic 3/pk



BOVIVET Vaginal Tape

For treatment of prolapsus of the uterus or vagina.

Cat. No	Description
152128	Vaginal tape for sheep, roll, 4 mm x 5 m



Vulva Suture Needles

For use with vaginal tape.

Cat. No	Description
152110	Gerlach suture needle, stainless steel, 19 cm
152120	Bühner suture needle, stainless steel, 30 cm



Plexi Speculum

For rectal and vaginal examinations.

Cat. No	Description
291000	Speculum plexi 20 x 240 mm with light
290979	Vaginal speculum for sheep 20 cm



Lamb Feeder

Specially for use in baby lambs or weak lambs.
Complete set with 60 ml syringe and feeding tube.

Cat. No	Description
247540	Lamb feeder



Nursing Bottle Lamb

Ideal for new born or very small lambs.

Cat. No	Description
247545	Nursing bottle for lamb with teat and valve



Feeding Buckets and Teats

10 l, the buckets are complete with valve, hanger and feeding teats.

Cat. No	Description
170823	Plastic bucket with 5 feeding teats
170822	Teats for 170823



Hoof Knives

Broad-bladed stainless, length of blades 80 mm. Standard quality.

Cat. No	Description
220052	Hoof knife standard right
220053	Hoof knife standard left
220054	Hoof knife standard double edge



Hoof Cutters

For sheep, 23 cm. Heavy, duty model.

Cat. No	Description
220224	Hoof cutter, toothed, for sheep
220225	Hoof cutter, for sheep



Sheep Shears

Standard sheep shears for wool cutting.


Cat. No	Description
247010	Sheep shears 29.5 cm, 6/pk



Tally Counter

An easy way to count your animals. 4 digits, counts 0-9999, chrome-plated.

Cat. No	Description
247135	Tally counter



Applicator and Castration Rings

For docking and tailing of lambs.

Cat. No	Description
200170	Applicator for castration rings
200180	Castration rings 100/pk, 50/pk



Trocar

The trocar ensures quick relief of pressure of rumen gas.

Cat. No	Description
190368	Trocar 4 x 120 mm with metal handle and cannula



INDEX

A	
AGRO examination gloves	26
Aluminium tube for sheath	100
Aluminium Wound Spray	52
Anti Kick Hobble	93
Antonius Farrowing Snare	115
Applicator and Castration Rings	122
Applicator for Suture Clips	77
Artificial Insemination Kit	98
Artificial Vagina Bull	98
Automatic Drenching Gun	47
Autopsy Knives	73
B	
BD Vacutainers and Needles	56
BHB milk test for Ketosis	7
Bleeding Knife	74
BLITZ Cartridges	79
BLITZ Slaughter Pistol	79
BLITZ Spare Parts	79
Blood transfusion Set	57
Blue-Spray	52
BoDia Quick	8
Bovine Uterine Metal Catheter	100
Brintoverilte	82
BTS Plus boar semen diluent	105
Bull Holder	93
Bull Ring Applicator	94
Bull Rings	94
C	
Caesarian Knife	73
Calf Drencher	65
Calf Fluid Feeder	66
Calfpullers	86
Calving Belt	85
Calving Rope	87
Carbon Steel Blades	73
Castration Forceps Pigs	117
Castrators for Bloodless Castration	78
Catheter Gel	100
Ceramic Magnets	69
Chlorhexidine Alcohol	82
Chlorhexidine Conc	82
Cleanser	81
Closing Caps for Syringes	40
CMT Test Kit	11
CMT Test Liquid	11
CMT Test Paddle	11
Coaxer	112
Colombus Teat Plugs in Box	89
Colostrum Densimeter	11
Colostrum Densimeter	66
Curved Scissors	76
Curved Snare Director	85
D	
Defiren Disinfection	81
Demineralized water	106
Demotec 90	58
Demotec 95	58

Demotec Easy Bloc	59
Demotec Futura Pad	59
Demotec Hoof Check	59
Demotec Hoof Discs	61
Demotec Wooden Blocks	59
Digital Thermometer DIGI-Vet SC	14
Digital Thermometer w/flexible tip, waterproof	4
Digital Thermometer, watertight	4
DIGI-TEMP Express Thermometer	4
Disposable Apron	30
Disposable Apron	31
Disposable Apron, non-sterile	31
Disposable Boots	31
Disposable Coat	30
Disposable Coveralls	29
Disposable Needle with Aluminium Hub	36
Disposable Needles	34
Disposable Needles	36
Disposable Obstetric Gown	29
Disposable Plastic Gown	30
Disposable Polyethylene Gown, Green	29
Disposable Scalpels	72
Disposable Shoe Covers	32
Disposable Staplers	77
Disposable Syringe	37
Disposable Syringes Luer Lock	38
Disposable Syringes with Fixed Needle	40
Disposable Syringes with Silicone O-ring	39
Drenching Guns	46
Drench-Matic Drenching Gun	47
Durasperm Plus boar semen diluent	106
E	
Ear Notcher	117
ECO-MATIC® Automatic Multi-dose Syringe	43
ECO-MATIC® Draw-off Caps	44
ECO-MATIC® TWIN – Automatic self-filling syringe with two separate barrels	43
ECO-MATIC® with Bottle Mount Universal Collars	43
Electric Dehorers	92
Embryotome and Accessories	84
Eye Hook	86
F	
FDA Marking Spray	48
Feeding Buckets and Teats	121
Feeding Tubes	118
G	
Gas Dehorers	92
Gel Granulate	54
Gel	54
Gel, Sterile	99
Glutavac Test	10
H	
Hand Disinfecting Gel	80
Handles for Obstetric Chains	86
Heating Lamp	114
Henke-Ject® TBC Syringe	45
Hoof Bandage	60
Hoof Cutter	62
Hoof Cutters	121

H

Hoof Glue	58
Hoof Hammer	63
Hoof Knives	121
Hoof Knives	61
Hoof Knives	61
Hoof Knives	61
Hoof Rasp	62
Hoof Rasps	60
Hoof Tar	63
Hoof Trimmer	62
Hoof-Dressing	60
Hoof-Gel	60
Hot-Shot Coaxer HS 2000	113
Hot-shot Coaxer Sabre-Six	113
Hot-Shot Dura Prod	113
Hot-Shot POWER-MITE	112

I

Infusion	56
Infusor Combi	57
Injection needles	34
Inoculating Sterile Needle	12
Insemination Catheter for Gilts	108
Insemination Catheter with Handle	107
Insemination Catheter	101
Insemination Catheters	107
Inseminator Universal	100
Intra-uterine insemination catheter	108

K

Kamar Heat Detector	102
Ketostix Ames	7

L

Laboratory Thermometer5,	104
Lamb Feeder	120
Lambing Instrument	119
Latex Examination Gloves - Powder Free	15
Latex Examination Gloves – Powdered	15
Latex Surgical Gloves - Powder Free	20
Latex Surgical Gloves – Powdered	21
LDA Toggle Sutures	77,78
Lighter Gas for Tail Docker	116
Liquid Catheter Gel	108
Liquid Hand Soap	80
Littmann Stethoscope	5
Lubricated Catheters	107

M

Magnet Introducer	69
Mannitol Salt Agar	10
Marking Spray	49
Marking Stick	48
Metal Detector	69
Metal Speculum for Vaginal Examinations	13
Metal Stand for Disposable Gloves	18
Milk Catheter	91
Milk Sucking Preventer	91
Milking Gloves - Powder Free	17
Milking Tube	90

Mobile Storage Box Boar Semen	105
MSU 1 Ultrasound Scanner for Pigs	110
MSU 2 Ultrasound Scanner for Pigs	111
Müller Hinton Agar	10
MULTI-MATIC® Syringes	44

N

Nitrile Clinic Gloves - Powder Free	17
Nitrogen Container Forceps	102
Nitrogen Container TR	103
Non-slip Disposable Boot Cover	32
Nose Clamp	93,94
Nursing Bottle Lamb	121
Nylon Hoof Bandage	63

O

0-Bite and Hartshorn Spray	118
Obstetric Apron	28
Obstetric Combi Suit	26
Obstetric Gloves Pigs	118
Obstetric Gown	26
Obstetric Hooks	87
Obstetric Jacket and Trousers	27
Obstetric Overalls	27
Obstetric Set	84
Obstetric Sleeve	28
Obstetrical Chains	86
Ointment Spray	52
Omniplast53,	117
Ovi-mark Marking Sprays	119

P

Panascope Stethoscopes	6
PiDia Light	9
Pig Holder	115
Plastic Buckets With Feeding Teats	66
Plastic Sucking Preventer	91
Plexi Speculum for sheep	120
Plexi Syringe	42
PLUS Plexi Revolver Syringe	42
Polansky Speculum for Vaginal Examinations	13
PORCIVET Boar Spray	109
PORCIVET Shoulder Protector	112
Post Mortem Kit	74
Post Mortem Tools	75
Practice Coats	32
Prolapse Pins	88

R

ReCovin ACT Paste	97
ReCovin Calcium Paste	96
ReCovin Dosing Gun	97
ReCovin Ketosis Paste	96
ReCovin Phosphate Paste	97
ReCovin Rumen Acidosis Paste	96
ReCovin Rumen Starter	97
Refractometers	12
Refractometers	66
Rondo Measuring Tape	13
ROUX-REVOLVER® Syringes	44
Ruminal Magnet	68
Ruminal Magnet, Extra Powerful	68

S

Safety Container for Sharp Objects	36
Saw Wire Handle	85
Saw Wire	85
Scalpel Handles	73
Scanner Gel	112
Schick Injector Razor	78
Semen Bags	106
Semen Bottles	106
Semen Container GT	103
Semen Container Voyageur	103
Semen filters	104
Sharpening Tools	62
Sheaths for Universal Inseminator	101
Sheep Restrainer	119
Sheep Shears	122
Shoulder Protector	28
Silicone Teat Insert	88
Skin Staple Removal Forceps	77
Sleeve Cover	18
Socorex Accessories	46
Socorex Automatic Self-filling Syringes	45
Soft Examination Gloves with Shoulder Protector	25
Soft Examination Gloves	25
Sonosleeve Examination Gloves	24
Sorting Paddle	114
Sorting Paddle	114
Sorting Plates	114
Spiral Insemination Catheters	108
Sprague Stethoscope	5
Sterile Gel	109
Sterile Scalpel Blades	72
Stomach Probang System Mark II	71
Stomach Pump	71
Stomach Tube	70
Storage Cabinet Boar Semen	105
Straight Scissors	76
Straw Cutter	101
Straw Tweezers	101
Styrko Kick Stopper	93
Super Pro Pig Tail Dockers	116
Super Sensitive Examination Gloves with Shoulder Protector	23
Super Sensitive Examination Gloves	23
Super Touch Examination Gloves	24
Super Touch Gloves with Shoulder Protector	24
Superior Marking Stick	48
Surgical Forceps	75
Surgical Forceps	75
Suture Clips	76
Swiss hoof knife	61
Syringe Mount Strips	90

T

Tail Docker	116
Tally Counter	118,122
Tar Bandage	63
Tattoo Figures	49
Tattoo Forceps, Single Row for Pigs	49
Tattoo Hammer	49
Teat bandage	53
Teat bandage	90
Teat Dilator	88

Teat Dilators	89
Teat Dipper	91
Teat Drainage Tubes	91
Teat Plugs Long	89
Teat Wipe	88
Tensoplast	53
Thawing Unit for Straws	102
Thermocult incubator	12
Thermometer Covers,	5
Thinklabs One	6
Thygesen Stomach Tube for Calves	70
Thygesen Stomach Tube for Cows	70
Tooth Cutting Forceps	117
Tooth Grinder	115
Trocar	71
Trocars	71,122

U

Ultrasound Scanner for Pigs	109
Umbilical Cord Clamps	118
Urine Stips	7
Uterine Support	119

V

Vaginal Plexi Speculum	13
Vaginal Speculum Polansky	99
Vaginal Speculum	99
Vaginal Suture Tape	87
Vaginal Tape	120
VALUELINE Cohesive Flex Bandage	50
VALUELINE Disposable Syringes	38
VALUELINE Examination Gel	54
VALUELINE Examination Gloves	25
VALUELINE Surgical Gloves – Powdered	22
Vet-Flex	50
Vet-Plast Elastic Adhesive Bandage	51
Vet-Vat Cotton	64
Vinyl Examination Gloves - Powder Free	16
Vinyl Examination Gloves – Powdered	16
Vulva Suture Needles	87,120

W

Wall Dispenser	80
Water Bath, Semen Diluent	104
Wax Teat Plugs	89
Wound Plast Spray with Tar	51
Wound Plast Spray	51

Z

Zinc Ointment 200 ml	52
----------------------	----

